

## Protocol

### Part A: Reverse transcription polymerase chain reaction (RT-PCR)

The use of *Tth*-DNA polymerase, which is also a reverse transcriptase in the presence of manganese and accepts both RNA and DNA as template, enables the whole reaction to be performed as a “one-step/one-tube” RT-PCR analysis (6). This experimental approach includes the amplification of fragments only to a maximum of 1 kb and a relatively high error rate for DNA polymerase as a result of the manganese ion concentration (7).

This one-step procedure minimizes the risk of cross-contamination of the sample.

Reverse transcription can be performed at 60 ° C minimizing problems from RNA secondary structure and high G/C content!

#### 1. One step RT - PCR (One step reaction eliminates the risk of cross contaminations associated with two step RT-PCR.)

“One-step” RT-PCR means reverse transcription and amplification is mediated by one enzyme in one tube.

The “one enzyme/one step” *Tth* DNA polymerase assay uses bicine buffers containing  $Mn^{2+}$  ions that are compatible with both RT and subsequent PCR (6, 11). **NOTE:** For RT-PCR amplification (reverse transcription and amplification in one step), it might be necessary to determine concentration of  $Mn^{2+}$  first by testing  $Mn^{2+}$  concentrations from 1- 4mM for each reaction.

#### Setting up “one-step” RT-PCR Reaction

Prepare two mastermixes 25µL each before starting.

**Table 1: Mix I for “one-step” RT-PCR**

Component	Volume	Final concentration
dNTP-Mix (10mM each)	1.5µL	300µM
Forward primer	Var.	450µM
Reverse primer	Var.	450µM
Template RNA	Var.	Up to 1µg (in steps of 1ng; 10ng; 100ng; 1µg)
Sterile water	Up to 25µL	
<b>Total</b>	<b>25µL</b>	

**Table 2: Mix II for “one-step” RT-PCR**

Component	Volume	Final concentration
5X RT-PCR buffer	10µL	1-time
Mn(OAc) <sub>2</sub>	5µL	2.5mM
<i>Tth</i> DNA Polymerase	0.5-1.0µL	2.5 - 5 units
Sterile water	Var.	
<b>Total</b>	<b>25µL</b>	

**NOTE:** Combine Mix I and Mix II on ice and gently vortex the final mixture in a PCR-tube. Place tube in Thermal cycler and run below mentioned (page 4) cycling programme immediately.

**Table 1: Possible 1-tube RT-PCR Cycle Protocol**

Step	Time	Temperature	Comments
<b>RT-Reaction</b>			
Initial Denaturation	1-3 minutes	95 °C	
RT-Reaction*	30 minutes	60° - 70 °C *	
<b>PCR-Reaction</b>			
<b>Initial Denaturation</b>	0.5 - 1 min.	95 °C	
<b>Cycling (1-10 step)</b>			
Denaturation	30-60 sec.	94 °C	
Annealing	30-60 sec.	50-70 °C	Approximately 5 °C to 8 °C below Tm of primers
Elongation	45-90 sec.	72-74 °C	
<b>Cycling (11-&gt;30 step) **</b>			
Denaturation	30 sec.	94 °C	
Annealing	30 sec.	50-70 °C	Approximately 5 °C to 8 °C below Tm of primers
Elongation	45 sec.	72-74 °C	
<b>Final elongation</b>	7 minutes	72-74 °C	

\*depends on your primers, 70 - 75 ° C is the optimal reaction temperature for *Tth* DNA polymerase. *Tth* DNA polymerase is resistant to prolonged incubations (20 min Half - Life time at 95 ° C ) at high temperatures (94 ° C) and can therefore be used for PCR.

\*\* Number of cycles ranges from 30 - 50 in the literature. If the template is limited, increased cycle numbers may result in nonspecific product yield.

## Part B: Two-step RT-PCR (optional - material not supplied)

In the two-step protocol the reversed transcription abilities of the *Tth* polymerase are used to perform a “standard” reverse transcription reaction for synthesis of cDNA. For this step the reverse primer is used together with the special RT-buffer system plus  $Mn^{2+}$ . The RT reaction is followed by a PCR reaction by adding PCR buffer,  $MgCl_2$  and forward primer.

As the error rate of *Tth* is increased at higher  $Mn^{2+}$  concentration (5), a two-step reaction is recommended if PCR products shall be used for subsequent cloning and/or sequencing experiments.

### 2. Two step RT - PCR:

Table 3: “two-step” RT-PCR reaction mixture - Components for RT-reaction

Component	Volume	Final concentration
10X Rev. Transcription buffer	2 $\mu$ L	1-time
Mn(OAc) <sub>2</sub>	2 $\mu$ L	0.9mM
dNTP-Mix (10mM each)	0.4 $\mu$ L	200 $\mu$ M
Reverse primer	Var.	450 $\mu$ M
Template RNA	Var.	up to 200ng
Tth-Polymerase	0.8 $\mu$ L	4 units
Sterile water	up to 20 $\mu$ L	
<b>Total</b>	<b>20<math>\mu</math>L</b>	For RT-reaction denature at 95 °C for 1-3 minutes, then incubate at 60-70 °C for 10-30 min.

### 2. Two step RT - PCR:

Table 4: “two-step” RT-PCR reaction mixture - PCR-reaction

Component	Volume	Final concentration
10X PCR reaction buffer	8 $\mu$ L	1-time
dNTP-Mix (10mM each)	0.4 $\mu$ L	200 $\mu$ M
Reverse primer	Var.	450 $\mu$ M
Forward primer	Var.	150nM
EGTA (7.5mM)*	10 $\mu$ L	0.75nM
Tth-Polymerase	0.8 $\mu$ L	4 units
Sterile water	up to 80 $\mu$ L	
<b>Total</b>	<b>80<math>\mu</math>L</b>	Gently vortex and add the 80 $\mu$ L PCR mastermix to the 20 $\mu$ L rev. transcriptase mixture after reverse transcriptase reaction. <b>Total PCR volume: 100<math>\mu</math>L</b> <b>Continue cycling immediately!!!</b>

\*\*\* For a 0.5M EGTA stock: dissolve 19.2g EGTA in 70mL deionized water, adjust pH8.0 with NaOH (10M). Add deionized water to 100mL final volume. Filter sterile with 0.22 $\mu$ m and store at room temperature

Table 5: Two-step PCR Protocol

Step	Time	Temperature	Comments
<b>Initial Denaturation</b>	1-2 minutes	95 °C	
<b>Cycling (1-10 step)</b>			
Denaturation	30-60 sec.	95 °C	
Annealing	30-60 sec.	50-70 °C	Approximately 5 °C to 8 °C below T <sub>m</sub> of primers
Elongation	45-90 sec.	72-74 °C	
<b>Cycling (11-&gt;30 step)**</b>			
Denaturation	30 sec.	94 °C	
Annealing	30 sec.	50-70 °C	Approximately 5 °C to 8 °C below T <sub>m</sub> of primers
Elongation	45 sec.	72-74 °C	
<b>Final elongation</b>	7 minutes	72-74 °C	

\*\* Number of cycles ranges from 30 - 50 in the literature. If the template is limited, increased cycle numbers may result in nonspecific product yield.

## Part C: “Standard PCR

*Tth*-polymerase can be used like *Taq*-polymerase for a standard PCR. As *Tth* is more unspecific for nucleotides than the other DNA-polymerases, *Tth* can be ideally used in the case of incorporation of modified nucleotides (digoxigenin, biotin or fluorescein). In all cases  $MgCl_2$ -concentration as well as the enzyme concentration in the PCR reaction has to be optimized. Typical  $MgCl_2$  concentrations are in the range of 1 - 6mM. The supplied standard concentration is 1.5mM  $MgCl_2$ . Optimal enzyme concentrations are in the range of 1.0 - 5.0 units. The standard concentration is 2.5 units per reaction.

### Setting up PCR Reaction mixtures

Prepare two mastermixes 50 $\mu$ L each before starting.

**Table 6: Mix I for standard PCR**

Component	Volume	Final concentration
dNTP-Mix (10mM each)	0.4 $\mu$ L	200 $\mu$ M
Forward primer	Var.	400 $\mu$ M
Reverse primer	Var.	400 $\mu$ M
Template DNA	Var.	up to 0.5 $\mu$ g
Sterile water	up to 50 $\mu$ L	
<b>Total</b>	<b>50<math>\mu</math>L</b>	

**Table 7: Mix I for standard PCR**

Component	Volume	Final concentration
10X PCR-buffer	10 $\mu$ L	1-time
<i>Tth</i> DNA Polymerase	0.5-1.0 $\mu$ L	2.5 - 5 units
Sterile water	up to 50 $\mu$ L	
<b>Total</b>	<b>50<math>\mu</math>L</b>	

**NOTE:** Combine Mix I and Mix II on ice and gently vortex the final mixture in a PCR-tube. Place tube in Thermal cycler and run below mentioned cycling programme immediately.

**Table 8: PCR Cycling Protocol**

Step	Time	Temperature	Comments
<b>PCR-Reaction</b>			
<b>Initial Denaturation</b>	1 - 2 min.	95 °C	
<b>Cycling (1-10 step)</b>			
Denaturation	30-60 sec.	94 °C	Approximately 5 °C to 8 °C below $T_m$ of primers
Annealing	30-60 sec.	50-70 °C	
Elongation	45-90 sec.	72-74 °C	
<b>Cycling (11-&gt;30 step) *</b>			
Denaturation	30 sec.	94 °C	Approximately 5 °C to 8 °C below $T_m$ of primers
Annealing	30 sec.	50-70 °C	
Elongation	45 sec.	72-74 °C	
<b>Final elongation</b>	7 minutes	72-74 °C	

\* Number of cycles ranges from 30 - 50 in the literature. If the template is limited, increased cycle numbers may result in nonspecific product yield.