---

fon:  
+49 (0)7357 - 91 63 77  
fax:  
+49 (0)7357 - 91 63 78  
eMail:  
info@genaxxon.com  
internet:  
www.genaxxon.com

## Chemically competent versus electro-competent cells

### Advantage of chemical competent cells

- Easy transformation of cells (only one pipetting step).
- Due to less pipetting steps contamination risk and error rate are reduced.
- Reaction time: max. 25 minutes (very often only 5 minutes instead of 20 minutes on ice are enough).
- No cross contamination by spoiled electroporation cuvettes.
- No additional equipment necessary (no Electroporator).
- No expensive electroporation cuvettes necessary.
- For "Heatshock" a simple waterbath or heating block can be used (even hot water from the tap can be used).
- A variable and, if necessary a big sample volume can be used. (up to 10 µl sample is no problem at all)
- High ionic strength does not influence the procedure. (ligation reaction contain high salt concentrations!!!)
- No total loss of sample as with electroporation (short-circuit) can occur.
- Transformation efficiency of  $> 1 \times 10^8$  is enough for most, even more complicated / demanding applications (No need for expensive high competent cells ( $10^9$ )).

### Advantage of electro-competent cells

- Electrocompetent cells are available with higher transformation efficiency.
- Transformation process faster (about 5 minutes, if everything works well).

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