# DNase I

## Deoxyribonuclease I [E.C. 3.1.21.1]

From bovine pancreas

Product	Cat#	Package size
DNase I, min. 3000 U/mg (Kunitz)	M3028.0010	10mg
DNase I, min. 3000 U/mg (Kunitz)	M3028.0050	50mg
DNase I, min. 3000 U/mg (Kunitz)	M3028.0500	500mg

### Product description

Deoxyribonuclease I (DNase I) from bovine pancreas is an endonuclease (glycoprotein), which preferentially cleaves the phosphodiester bond in DNA behind pyrimidine nucleotides. This results in a polynucleotide with a 5'-phosphate and a free OH-group in position 3'. DNase I cleaves single-stranded and double-stranded DNA as well as chromatin. The specificity of the enzyme reaction (single-strand-“Nicks” versus double-stranded breaks) is determined by the ions available. In the presence of Mg<sup>2+</sup> single-strand nicks are generated and in the presence of Mn<sup>2+</sup> double-strand breaks. The pH-optimum of DNase I is 7.8 and it is activated by divalent cations. Maximum activation requires the presence of Mg<sup>2+</sup> and additional Ca<sup>2+</sup>. Calcium ions (5mM) protect DNase I from proteolytic digest. Inhibition is achieved by citrate, if activation is done by magnesium, but not if manganese has been the activator. Besides, it is inhibited by chelators such as EDTA and SDS or β-mercaptoethanol.

DNase I is used, among other things, to avoid unwanted cell clumping during tissue disaggregation, which makes dispersion of individual cells difficult. Since tissue disaggregation and subsequent isolation of single cells is always accompanied by rupture and lysis of some cells, DNA is released from these cells into the surrounding medium. The released DNA in the medium is thought to be responsible for undesirable cell clumping during tissue dissociation procedures. This clumping effect can be avoided by adding deoxyribonuclease (DNase) to the dissociation medium and has been described for a variety of different cell and tissue types.

The enzyme can also be used in molecular biology techniques like digestion of DNA, in the RNA purification (ref. 2 Suppl. 1 pp. 4.1.4) or generating “random nicks” for “nick translation” (ref. 2 Suppl. 9 pp. 3.5.4-6) or “footprint”-assays (ref. 2 Suppl. 7 chapter 12.4) or investigation on chromatin (ref. 2 Suppl. 48 chapter 21.4.1).

**NOTE:** This product is NOT recommended for isolation of RNA without prior chemical treatment to get rid of remaining RNase activity as this enzyme preparation is not guaranteed to be free of RNase activity.

### Preparation of stock solution

DNase I is used to digest DNA in the isolation and purification of RNA. Dissolve DNase at 2mg/mL in 0.15M sodium chloride or reaction buffer (e.g., 50mM Tris/HCl pH7.5, 10mM MgCl<sub>2</sub> (single strand nicks) or 10mM MnCl<sub>2</sub> (double strand breaks) respectively, 50µg/mL BSA).

### Storage and stability

For stability reasons the concentration of DNase I containing solutions should be at least 1mg/mL. Do not vortex.

This solution is stable for more than one year at -20°C (ref. 2 Suppl. 8 page 3.12.5).

The lyophilised form is stable for 2-5 years if kept at +2°C to +8°C.

If a solution is protease-free, DNase I will not lose significant activity at pH5-7 and 62°C for 5 hours.

**Heat-inactivation of DNase I:** Heat to 99°C for 10 minutes.

### Unit definition

One unit is defined as the amount of enzyme required to completely degrade 1µg of lambda DNA in 10 minutes at 37°C in 50µL of a buffer containing 40mM Tris/HCl (pH8.0), 10mM NaCl, 6mM MgCl<sub>2</sub> and 10mM CaCl<sub>2</sub>. Under these conditions one unit of DNase I activity is comparable to 1 Kunitz unit.

**Alternative definition:** One unit is defined as the amount of enzyme which causes an increase of absorbance at 260nm of 0.001 per minute at 25 °C and pH5.0 based on the method of Kunitz.

### Preparation of RNase-free DNase I (according to ref. 2, page 3.12.6 Supplement 8):

Since DNase I may contain trace amounts of RNases, it is recommended that DNase I is treated with the following protocol prior to use.

- Dissolve 1mg/mL DNase I in 0.1M iodoacetic acid plus 0.15M sodium acetate at a final pH of 5.3.
- Heat to 55°C for 40 minutes, then cool down the solution.
- Add 1M CaCl<sub>2</sub> to a final concentration of 5mM.
- Store in small aliquots at -20°C.

### Material needed but not supplied:

- 10X Reaction Buffer: 400mM Tris/HCl (pH8.0), 100mM MgSO<sub>4</sub> and 10mM CaCl<sub>2</sub>.
- Storage Buffer: 20mM Tris/HCl (pH7.5), 50% glycerol (v/v), 1mM MgCl<sub>2</sub>. **NOTE:** conc. of DNase I: not below 1mg/mL!
- Ca<sup>2+</sup> and Mg<sup>2+</sup> or Mn<sup>2+</sup>
- Stop solution: 20mM EGTA (pH8.0)

### Associated activities

To test RNase activity, 50ng of [<sup>3</sup>H]RNA is incubated with 5 units DNase I in transcription buffer for 1 hour at 37°C, and the release of radiolabeled nucleotides is monitored by scintillation counting of TCA-soluble material. For passing the test the maximum count has to be below 3%.

### Inhibitors

EGTA; EDTA; DTT; salt concentration >100mM

### Specifications

<b>Source:</b>	bovine pancreas
<b>Molecular weight:</b>	approx. 31 kDa
<b>CAS number:</b>	[9003-98-9]
<b>Activity:</b>	min. 3000 U/mg (Kunitz)
<b>Solubility:</b>	Dissolves readily at 5mg/mL in 150mM NaCl to give a clear, colourless solution.
<b>RNase activity:</b>	not guaranteed free of RNase activity
<b>Storage:</b>	-20°C

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