Proteinase K

>30 units/mg (MB grade)

incl. Application protocols

inkl. Anwendungsprotokolle

Cat#: M3036
Cat#: M3037

Contact & Technical support

Tel.: 0731 3608 123

Fax: 0731 3608 962

E-Mail: info@genaxxon.com

Description:	Proteinase K is a non-specific serine protease having a very high specific activity. It has been used for isolation of mRNA, high molecular weight DNA and to inactivate other enzymatic activities. Proteinase K is active with or without the presence of SDS and EDTA.
Source:	<i>Parengyodontium album</i> (<i>Tritirachium album</i>).
Host:	<i>Komagataella phaffii</i> (<i>Pichia pastoris</i>).
Activity:	>30 units/mg lyophilizate (Molecular Biology (MB) grade)
Specific activity:	>40 units/mg protein (Molecular Biology (MB) grade)
Form:	lyophilised powder (M3036) frozen solution (TE, pH8.0) of 20mg/mL (>600 units/mL) (M3037)
Unit definition:	One unit is defined as the amount of enzyme that liberates Folin-positive amino acids and peptides, corresponding to 1 pmol tyrosine under assay conditions in 1 minute using haemoglobin as substrate (1).
Stability and storage:	Proteinase K is stable at RT for short periods of time. For long term storage, we recommend -20°C.

Preparation of Proteinase K solution: The standard concentration of a Proteinase K solution is 20mg/mL. For preparing this solution, 20mg of Proteinase K is dissolved in 1mL pure (de-mineralised, DNase and proteinase-free) water. The solution is stable at -20°C for 2 years but can also be stored at +2°C to +8°C for several months.

QUALITY CONTROL

16-hour incubation:	a 50µL reaction solution containing 1µg of lambda-DNA and 1.8 units enzyme incubated for 16 hours at 37°C resulted in the same DNA band pattern after gel electrophoresis as compared to the pattern produced without enzyme.
Exonuclease activity:	Incubation of 6 units of the enzyme for 4 hours at 37°C in 50µL assay buffer with 1g sonicated 3H DNA (3×10^5 cpm/µg) released less than 0.2% of radioactivity.
Endonuclease activity:	Incubation of 1.8 units of enzyme with 1µg PhiX174 RFI DNA in 50µL assay buffer for 4 hours at 37°C gave less than 1.5% conversion of RFI.
RNase contamination:	Incubation of 6 units of enzyme with 1 µg MS2 RNA in 50µL assay buffer for 4 hours at 37°C resulted in the same RNA band pattern after gel electrophoresis as compared to the pattern produced without the enzyme.

(1) Anson M.M (1939) J. Gen. Physiol., 22, 79.

Product	Cat#	Package size
Proteinase K (from <i>Tritirachium album</i>): >40 units/mg protein	M3036.0100	100mg
Proteinase K (from <i>Tritirachium album</i>): >40 units/mg protein	M3036.0200	200mg
Proteinase K (from <i>Tritirachium album</i>): >40 units/mg protein	M3036.0500	500mg
Proteinase K solution (20mg/mL): >600 units/mL	M3037.0001	1mL
Proteinase K solution (20mg/mL): >600 units/mL	M3037.0005	5 x 1mL

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Proteinase K (English description)

lyophilised (>30 U/mg lyophilizate - MB grade)
isolated from *Pichia pastoris*

Proteinase K belongs to the family of subtilisin-like serine proteases with very strong proteolytic activity. The enzyme is secreted from *Tritirachium album* cultures, whereas the "K" stands for the Kreatin (German for Creatine) the main source of carbon and nitrogen source of *Tritirachium album*.

Proteinase K is an unspecific protease that shows endo- and exoproteolytic activity. Activated by calcium (1 - 5mM), the enzyme digests proteins preferentially after hydrophobic amino acids (aliphatic, aromatic and other hydrophobic amino acids). Proteins will be completely digested, if the incubation time is long and the protease concentration high enough.

The enzyme has two binding sites for Ca²⁺ ions which do not really participate in the catalysis but do contribute to the structural stability of the protein. Upon removal of the calcium ions, the stability of the enzyme is reduced, but the proteolytic activity remains. Proteinase K has two binding sites for Ca²⁺, which are located close to the active center, but are not directly involved in the catalytic mechanism. Removal of the Ca²⁺-ions reduces the catalytic activity of proteinase K by 80 %. The residual activity is sufficient to digest proteins, which usually contaminate nucleic acid preparations. Therefore, the digest with proteinase K for the purification of nucleic acids is performed in the presence of EDTA (inhibition of magnesium-dependent enzymes).

If the presence of Ca²⁺ is required, Ca²⁺ is added up to a concentration of 1mM and is removed by the addition of EGTA (pH8,0; final conc. 2mM) later on. The pH-optimum is at pH8, but the enzyme is active over a wide pH-range (pH4,3 - 12). An elevation of the reaction temperature from 37°C to 50 - 60°C may increase the activity several times, like the addition of 0,5 - 1% SDS. Temperatures above 65°C, trichloroacetic acid or the serine protease-inhibitors AEBSF, PMSF or DFP will inhibit the activity. Proteinase K will not be inhibited by EDTA, urea (1 - 4M), SDS, citrate, iodoacetic acid or, by other serine protease inhibitors like TLCK and TPCK.

Proteinase K can be inactivated by incubation at a temperature of 95°C. Make sure that inactivation temperature is not below 95°C and inactivation time not shorter than 10 minutes.

Proteinase K NGS Grade is developed for most demanding applications. Additional purification technology results in its significantly increased solubility (2.5-fold), increased specific activity, and remarkable purity with DNA content ≤0.1pg/mg, compared to our Molecular Biology Grade product.

Proteinase K is used for the destruction of proteins in cell lysates (tissue, cell culture cells) and for the release of nucleic acids since it inactivates DNases and RNases very effectively.

Following some application examples:

Based on 20mg/mL concentrations. Proteinase K dissolved in pure water or 50mM Tris, pH8.0.

Application Protocols:

- Isolation of genomic DNA from mammal cells
(3h, 50 °C) 10mM Tris/HCl (pH8.0), 100mM EDTA (pH8.0), 50mM NaCl, 0.5%SDS, 20µg/mL DNase free RNase, 100µg/mL proteinase K.
- Isolation of genomic DNA from mice tails
(overnight, 55 °C) 20mM Tris/HCl (pH8.0), 5mM EDTA (pH8.0), 400mM NaCl, 1% SDS, 400µg/mL proteinase K.
- Rapid isolation of genomic DNA as a PCR-template from cell cultures
(1h, 37 °C) 67mM Tris/HCl (pH8.8), 16.6mM ammonium sulfate, 5mM β-mercaptoethanol, 6.7mM MgCl₂, 6.7µM EDTA (pH8.0), 1.7µM SDS, 50µg/mL proteinase K.
- Preparation of DNA for PFGE, lysis in agarose block
(approx. 20h, 50 °C). 100mM EDTA (pH8.0), 1 mM Tris/HCl (pH7.6), 20mM NaCl, 1% sarcosyl, 100µg/mL proteinase K
- Purification of mRNA prior to cDNA synthesis
(often leads to an increase in cDNA yield)
(1-2h, 37 °C) 100mM Tris/HCl (pH7.5), 10mM EDTA (pH8.0), 50mM NaCl, 0.1% SDS, 5µg/mL proteinase K.
- Purification of PCR-reactions prior to cloning
(overnight, 55 °C) PCR-reaction + 10mM Tris/HCl (pH8.0), 5mM EDTA (pH8.0), 400mM NaCl, 1% SDS, 400µg/mL proteinase K.
- Purification of plasmid-DNA prior to in vitro transcription
(1h, 37 °C) Plasmid-DNA + 21mM Tris/HCl (pH8.0), 5mM EDTA (pH8.0), 50mM NaCl, 0.5% SDS, 100µg/mL proteinase K.
- RNase inactivation in ribonuclease-protection-assays
(30 min, 37 °C) 30µL RNA/RNA-hybridisation mix + 300µL RNase digestion mix + 0.6% SDS, 300µg/mL proteinase K.
- Inactivation of DNase in transcription-run-on-assays
(30 min, 42 °C) Suspension of nuclei in 36mM Tris/HCl (pH7.4), 17mM MgCl₂, 0.7mM CaCl₂, 170mM NaCl, 13mM EDTA (pH8.0), 0.5% SDS with DNase, 100µg/mL proteinase K.
- Denaturation of alkaline phosphatase when capping vector DNA
(30 min, 56 °C) 1 x dephosphorylation buffer + 0.5% SDS, 5mM EDTA (pH8.0), 100µg/mL proteinase K.

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