Genaxxon BioScience

N-terminal Protein Sequencing Service

fon:
 +49 (0)7357 - 91 63 77
 fax:
 +49 (0)7357 - 91 63 78
 eMail:
 info@genaxxon.com
 internet:
 www.genaxxon.com

Product	Cat#	Package size
Automated Edman Degradation of Protein/Peptide	P2047.0000	Set-up
Subsequent Degradation steps	P2048.0000	step
Chemical or enzymatic cleavage of protein sample in solution	P2049.0000	Set-up
Proteolytic cleavage of protein/peptid sample	P2050.0000	Set-up
Fractionation of degradation peptides by HPLC	P2051.0000	Set-up
Insufficient sample amount (set-up fee)	P2053.0000	Set-up

Product description

For being able to N-terminal sequence customers proteins or peptides Genaxxon needs lyophilized samples, in solution or blotted on PVDF membrane (blotting can be provided by Genaxxon as fee-for-service). Proteins and peptides can be sequenced five to ca. 40 steps depending on sample amount, amino acid sequence and customer's requirements.

As result customer will receive a report of the N-terminal sequence analysis with the detected amino acid sequence and all HPLC chromatograms as a PDF file by email or as printed version.

Application area for protein sequencing

- * Determination the N-terminal amino acids of a protein
- * Determination the N-terminal amino acids of a peptide
- * Check of the correct translation of a recombinant protein
- * Check of the sequence and purity of a recombinant protein
- * Check of the sequence and purity of a synthetic peptide
- * Determination of complete protein sequence in combination with additional methods
- * Quality control

How to order

Please contact us at first. We can give you some helpful information and clarify open question points. Then ship sample together with additional information (as clarified before) to us.

Prices N-terminal sequencing

Service - N-terminal protein and peptide sequencing	Net Preis (without VAT)
Basic service (Set-up with first 5 amino acid sequencing steps)	367.50 Euro
Set-up of analysis / minimal price (e.g. if no sequence was detected due to N-terminal blockage)	273.00 Euro
each additional sequencing step	52.50 Euro
Discount for at least 10 protein/peptide samples	10 %
Extra charge for Express Service (The analysis will be performed within 3-7 days after arrival of sample. Please contact us at first!)	50% extra charge to the basic service (incl. 5 steps); 25% extra charge for each additional sequencing step

Sample Requirements

Amount		Min. 20pmol. Smaller amounts are only accepted after detailed discussion!
Sample type		<ol style="list-style-type: none"> 1. lyophilized protein/peptide samples 2. liquid protein/peptide sample in 10-100µL volatile solvent like Milli-Q water, isopropanol, acetonitrile. The sample have to be shipped in frozen state! 3. Protein samples blotted on PVDF membrane (PVDF membrane size: max. 3x6 mm, smaller and more concentrated sample are preferred. The protein band can be stained by Ponceau S, Amino black, Sulforhodamine or Coomassie Brilliant Blue.
Purity		The sample should be as pure as possible (at least 75-80% purity) and contain only the protein or peptide. Free amino acid, primary amines, SDS, salt, buffer or other contaminants should be as low as possible since the Edman chemistry can be negatively affected.
Cysteine modifikation		Cysteine without special modification can not be detected by N-terminal sequencing. Therefore the sample has to be modified for detection of cysteine. Below you can find a modification protocol.
N-terminal blockage		Proteins and peptides which are N-terminally blocked do not have a free N-terminal amino group. Therefore these proteins can not be sequenced. More than 50% of all eukaryont proteins are blocked. On request we can perform deblocking procedures but we need a significant higher amount of protein and the deblocking does not work always, because mostly the type of blockage is unknown.
Glycosilation and other modifications		Edman sequencing steps without the detection of a PTH amino acid, reduced peak intensity or altered retention times can be caused by glycosilation, phosphorylation or other modifications. By mass spectrometry the modification can be analyzed.

Details for sample shipping to Genaxxon

- * Support Genaxxon with information about the sample (e.g. sample name, sample type, MW, protein amount, staining method) and the requested number of N-terminal sequencing steps. Sign details.
- * Send us an email or call us before shipping the samples
- * The samples should be in a reaction tube sealed by Parafilm. Protein samples on PVDF membrane can either be shipped as small and dry membrane slides in reaction tubes or as a dried PVDF membrane in plastic foil with enclosed description of the protein bands which should be analysed. We cut the marked protein band for analysis.
- * Liquid samples should be sent in frozen state for precaution.
- * Lyophilized sample and samples blotted on PVDF membrane can typically shipped at room temperature.
- * Please, send your samples in a padded envelope or in a box together with the Edman protein sequencing details.

Our Recommendation for Protein Blotting

Semidry PVDF Blotting Protocol

We recommend the following protocol for semidry blotting on PVDF membrane:

- PVDF membrane: Immobilon P membrane or comparable
- Blot buffer: 50 mM sodium borate, pH 9.0 / 20% methanol (HPLC quality)
(0.1% SDS can be added to blotting buffer if protein above 40kDa should be sequenced)
- Blotting condition: 1mA/cm1 PVDF membrane for 2-3 hours at 4°C

Protein staining of PVDF membranes

Ponceau Red Staining

Staining solution (0.5% Ponceau S, 1% acetic acid in Milli-Q water):

- 0.25g Ponceau S
- 0.5mL acetic acid
- ad 50mL Milli-Q water

Procedure:

1. Wash the PVDF blot membrane 2 x 3 minutes with plenty Milli-Q water.
2. Stain the PVDF membrane with Ponceau S staining solution for 1-3 minutes.
3. Destain the PVDF blot membrane under visual control with Milli-Q water until protein bands are well visible.
4. Dry the PVDF membrane.

Sulforhodamine Staining

Staining solution (0.005% sulforhodamine, 0.2% acetic acid, 30% methanol in Milli-Q water):

- 150mL methanol
- 1mL acetic acid
- 25mg sulforhodamine
- ad 500mL Milli-Q water

Procedure:

1. Wash the PVDF blot membrane 2 x 10 minutes with plenty Milli-Q water.
2. Dry the PVDF membrane at room temperature!
3. Stain the PVDF membrane in Sulforhodamine staining solution for 1-2 minutes.
4. Wash the PVDF membrane with Milli-Q water shortly.
5. Dry the PVDF membrane.

Coomassie Blue R250 Staining

Staining solution (0.1% CBB R250, 10% acetic acid, 40% methanol in Milli-Q water):

- 0.1g CBB R250
- 40mL methanol
- 10mL acetic acid
- ad 100mL in Milli-Q water

Destaining solution (10% acetic acid, 40% methanol in Milli-Q water):

- 40mL methanol
- 10mL acetic acid
- ad 100mL Milli-Q water

Procedure:

1. Stain the PVDF blot membrane for 5 minutes in CBB R250 staining solution
2. Destain the PVDF blot membrane for 3 x 5 minutes with destaining solution under visual control until protein bands are well visible.
3. Dry the PVDF membrane.

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