

RNA Isolation Protocol

1. **Place** the fragmented biological material in a 2mL tube.
Add 600µL RLys Buffer and vortex for 60 seconds
2. **Centrifuge** for 120 seconds at $\geq 12000 \times g$ (preferably at $15000 \times g$).
3. **Transfer** supernatant into an **RNase-free** 1.5 or 2.0mL reaction tube and add **600µL 70% ethanol** to the transferred supernatant.
 Mix by pipetting or vortexing.
NOTE: For homogenization using bead-beating tubes: carefully pipet the appropriate volume of the supernatant by placing a 200µL pipette tip (N.B.: a 1mL tip may be clogged by the beads) into the filling. Tissue remains should either lie on one side of the tube or at the bottom.
4. **Transfer** up to **700µL of the obtained mixture** into an **RNA Purification Column** placed in a collection tube. Centrifuge for 15 seconds at $\geq 12000 \times g$.
 Discard flow-through and re-use the column, together with the collection tube.
5. **Transfer** the **remaining mixture** into the same purification column and centrifuge for 15 seconds at $\geq 12000 \times g$.
 Discard flow-through and place the mini column in a new 2mL collection tube.

6. Optional - On-column DNase treatment

- a. Prewash mini column with 500µL **RW2 Buffer** and centrifuge for 60 sec. at $\geq 12000 \times g$.
 Discard flow-through and re-use the collection tube.
 - b. For each isolation mix 90µL **DNase I Reaction Buffer** and 10µL reconstituted **DNase I (not included in this kit)**. Mix by inverting the tube.
 - c. Apply 95µL of the above mixture onto the centre of the RNA Purification Column.
 Incubate for 5 minutes at room temperature.
 - d. Add **600µL RW1 Buffer** and centrifuge for 15 sec. at $\geq 12000 \times g$.
 Discard flow-through and re-use the collection tube.
Proceed with step 8.
-

7. **Omit** this step in the case of DNase I treatment / step 6)
Add 700µL RW1 Buffer and centrifuge for 15 sec. at $\geq 12000 \times g$.
 Discard flow-through and reuse the collection tube.
NOTE: Do proceed with step 8 directly without step 7 if the optional DNase I treatment was done, as described in step 6.
8. **Add 500µL RW2 Buffer** and centrifuge for 15 sec. at $\geq 12000 \times g$.
Repeat this step.
9. Centrifuge for 90 sec. at $\geq 12000 \times g$ (preferably at $15000 \times g$). Discard the collection tube and the flow-through and carefully transfer the Purification Column to a RNase-free 1.5mL reaction tube.
NOTE: **RW2 Buffer** contains alcohol, which may interfere with some enzymatic reactions and decrease the elution efficiency. It is therefore crucial to remove the alcohol completely from the mini column before elution.
10. **Add 50-100µL elution buffer REB** precisely onto the centre of the purification column membrane. Centrifuge for 60 sec. at $\geq 12000 \times g$ to elute purified RNA.
 The isolated RNA is ready for use in downstream applications or for storage at -80°C .
NOTE: Other buffer volumes in the range of 30-50µL may be used. For instructions, see page 8 (RNA elution) (Recommendations and important notes).

RNA Clean-up Protocol (for already isolated RNA)

1. Adjust RNA sample volume to **100µL** by adding RNase-free water to a 1.5 or 2.0mL RNase-free reaction tube.
2. **Add 300µL RLys Buffer.**
3. **Add 300µL 96-100% ethanol.** Mix well by pipetting or vortexing.
4. Transfer the mixture into a **RNA Purification Column** placed in a collection tube.
 Centrifuge for 15 sec. at $\geq 12000 \times g$.
 Discard the flow-through and re-use the collection tube.
 Proceed to **step 7** of the Isolation Protocol (page 7)