Genaxxon Bioscience User Guide for Polycarbonate Capillaries

Manual for using the Genaxxon
polycarbonate capillaries
in Roche LightCycler® machines

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Product	Cat#	Package size
Polycarbonate Capillaries for use in capillary PCR	I2250.0960	960 capillaries
Starter Kit with polycarbonate capillaries	I2250.0961	192 capillaries + 1 carousel 1.5
Starter Kit with polycarbonate capillaries	I2250.0962	192 capillaries + 1 carousel 2.0
PEEK-Carousel for use in LightCycler™ 1.2/1.5	I2250.0001	1 carousel 1.2/1.5
PEEK-Carousel for use in LightCycler™ 2.0	I2250.0002	1 carousel 2.0

Hotline: +49 89 70009551 or info@genaxxon.com

Notes on Warranties and Disclaimer

Genaxxon is dedicated to your success and every batch of this product is tested with an extensive routine procedure to make sure that it meets all your needs. However, it has neither been developed nor tested for a specific application.

This product is for research use only. For *in vitro* use only

Genaxxon's liability with respect to any product is limited to the replacement of the product. No other warranties are provided by Genaxxon. Genaxxon is not liable to any direct, indirect, incidental or consequential damage arising out of or in connection with the use of any of Genaxxon products.

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Kit Contents & Storage Information

I2250.0961 contains 1 carousel for use in the LightCycler® Version 1.2/1.5.
2 x 96 capillaries and caps and 1 transfer-pin for handling the caps.

I2250.0962 contains 1 carousel for use in the LightCycler® Version 2.0.
2 x 96 capillaries and caps and 1 transfer-pin for handling the caps.

Capillaries can be used for either LightCycler® Version, respective either carousel.

Capillaries are certified DNase, RNase and Protease free.

Capillaries, carousels and transfer-pin can be stored for long term at Room Temperature without loss of their usability.

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Description

Sample preparation / PCR Mixture preparation

Samples for PCR can be handled and prepared according to the existing protocol at customer site. No changes have to be made using polycarbonate capillaries instead the original glass capillaries from Hofmann La-Roche.

No changes have to be made preparing the PCR Mix. All PCR-Mastermixes that work in glass capillaries will also work in polycarbonate capillaries.

It is possible to use our polycarbonate capillaries for PCR volumes from 10 µl to 50 µl without changes in the PCR programme (for details see point PCR programme).

Handling capillaries and caps

For better/easier handling of the capillaries and caps we provide a "Transfer-Pin" for taking the caps out of the rack-box. The Transfer-Pin is part of our "Starter-Kit". If you are used to work without a "Transfer-Pin", you can also work without it, but use protective gloves to protect you PCR-reaction from contamination, by for example DNAses.

For sealing the capillaries with the caps, please use the transfer pin, to transfer caps from the rack-box to the polycarbonate capillary. Ensure that each capillary is closed tightly by pressing the caps tightly into the capillaries.

To ensure that the capillaries fit in the cavities of the carousel, the capillaries have to be pressed into these cavities. In opposite to the glass capillaries from Roche, it will not harm the capillaries nor break them.

The Carousels

Genaxxon offers 2 types of carousels that fit to LightCycler® Version 1.2/1.5 or to Version 2.0 respectively. The carousel type 1.2/1.5 will only fit into the LightCycler® Version 1.2/1.5. It will not be possible to use it in the LightCycler® Version 2.0. As this carousel weights about 20 g more than the 2.0 Version it is not recommended to use it in the centrifuge together with the centrifuge rotor for the LightCycler® carousel Version 2.0.

NOTE: Genaxxon is not liable for any damage that occurs for the wrong usage of the carousel in the centrifuge.

The LightCycler® carousel 2.0 will fit in the LightCycler® 1.2/1.5 as well as the LightCycler® 2.0. As this carousel is about 20 g lighter than the 1.2/1.5 type it can only be used in the centrifuge with the appropriate rotor for the LightCycler® 2.0 carousel!

NOTE: Genaxxon is not liable for any damage that occurs for the wrong usage of the carousel in the centrifuge.

Please inform Genaxxon about the LightCycler® version available in your laboratory.

Description (continued)

Using LightCycler® centrifuge adaptors

PCR master mixes and additional components can be easily pipeted into the upper part of the capillaries (exactly as with the glass capillaries from Roche). To get the liquid to the bottom of the capillaries, the capillaries can be put into a 1.5 ml reaction tube (cut caps or close them) and then be transferred to a table top centrifuge (very short centrifugation (1-2 sec.) is enough to get the liquid down). An even easier possibility is to “shake” the liquid down with a short jerk.

In either cases it is not necessary to use the centrifuge adaptors from Roche. If you want to use the adaptors, this will be possible without changes in handling compared to the glass capillaries from Roche. The Genaxxon capillaries however do not fit completely in the Roche centrifuge adaptors (the capillary head will stick out). For that reason it will not be possible to put the rotor lid on top of the centrifuge rotor of mini centrifuges. As a very short centrifugation is enough to transfer the liquid from the top of the capillaries to the bottom is not necessary to put the rotor lid on the rotor.

Performing PCR / PCR programme

The Genaxxon polycarbonate capillaries have a bigger diameter and are more thick walled than the 20 µl glass capillaries from Roche. For this reason and to guarantee effective heat transfer from the outside of the capillaries to the middle of the PCR reaction, it might be necessary to increase the “hold-times” at Denaturing and Annealing.

For this reason we recommend strongly to program your LightCycler® with at least 10 seconds at 95°C for the denaturing step and 10 seconds for the annealing step. Too short denaturing and annealing times will lead to a complete failure of your PCR, or at least to higher Ct/Cp values. According to the specific PCR it might even be necessary to use a denaturing time of 15 seconds.

If you are already using a programme with both times higher than 10 seconds, it should not be necessary to make changes in the program. Nevertheless you can try to increase the existing denaturing and/or annealing times by additional 5 seconds. If this will not lead to better results it makes no sense to increase denaturing and/or annealing times further.

As the easiness of denaturing DNA depends on the GC-content of your DNA, it might be better to increase the denaturing temperature from 95°C to 97°C (for not longer than 10 sec.). Longer times will not improve results.

NOTE: Longer denaturing times at 97°C will rarely give better results.

Some examples of temperature protocols are shown on the next page.

Example (95°C Denaturing)

	Glass	Polycarb.	Glass	Polycarb.	Glass	Polycarb.
Denaturing	0 sec.	10 sec.	5 sec	10-12 sec	> 7 sec	> 10 sec
Annealing	0 sec.	10 sec	5 sec	10-12 sec	> 7 sec	> 10 sec
Extension	10 sec.	10 sec.	5 sec	5-8 sec	> 7 sec	> 7 sec

Example (97°C Denaturing)

	Glass	Polycarb.	Glass	Polycarb.	Glass	Polycarb.
Denaturing	0 sec.	5 sec.	5 sec	7-10 sec	> 7 sec	< 10 sec
Annealing	0 sec.	10 sec	5 sec	10 sec	> 7 sec	> 10 sec
Extension	10 sec.	10 sec.	5 sec	5-8 sec	> 7 sec	> 7 sec

Fluorescence Detection

All standard dyes (probes) can be detected while using the Genaxxon polycarbonate capillaries.

It might be a problem to detect ROX in the PCR reaction. If you use ROX and you get no PCR signal at all, perform a standard gel electrophoresis. If you can see a DNA-band, but you get no Real-Time-PCR signal it will not be possible to improve the PCR by changing temperatures or hold-times.

Melting Curve Analysis

If the melting curves show higher values or are spread over a wider range while analysing a dilution series, it is recommended to decrease the slope of the heating step (less than 0.1 increments).

Sample recovery after PCR

If you want to analyse the PCR reaction on a gel it is the easiest way to cut the capillaries with scissors and spin the solution into a clean reaction tube from which you can pipet the PCR reaction to a gel or what ever purpose you need the PCR reaction for.

It will be also possible to remove the caps and put the capillaries upside down into a reaction tube spinning the tube with the capillary.