

Ribonuclease A - DNase-free

DNase-free RNase A (E.C. 3.1.27.5)

DNase-free Ribonuclease A (E.C. 3.1.27.5) from bovine pancreas (salt free, freeze dried)

Product	Cat#	Package size
RNase A – Dnase-free, min. 80 U/mg (Kunitz)	S5218.0050	50mg
RNase A – Dnase-free, min. 80 U/mg (Kunitz)	S5218.0250	250mg
RNase A – Dnase-free, min. 80 U/mg (Kunitz)	S5218.0500	500mg

Product description

Ribonuclease A- DNase-free is a DNase and Protease-free enzyme that can be used directly for RNA digests.

Ribonuclease A (RNase A) is an endoribonuclease, that specifically cleaves single-stranded RNA 3' to pyrimidine residues (cytosine, uracil). Thereby, it generates pyrimidine-3'-phosphate or oligonucleotides with terminal pyrimidine-3'-phosphates. For example the sequence of pG-pG-pC-pA-pG will be cleaved to give pG-pG-pCp and A-pG. The highest activity is exhibited with single stranded RNA. The pH-optimum is in the range of 7.0 to 7.5.

RNase A is a single chain polypeptide containing 4 disulfide bridges. In contrast to RNase B, it is not a glycoprotein.

RNase is used for the purification of RNA-free DNA, for the removal of non-hybridised regions of RNA : DNA hybrids or as a molecular weight marker. The enzyme is inhibited by diethyl pyrocarbonate (DEPC), guanidinium salts (4M GuaSCN), β -mercaptoethanol, heavy metals, vanadyl-ribonucleoside-complexes, RNase-inhibitor from human placenta and competitively by DNA while the effect of denatured DNA is higher than by native nucleic acids, or by alkylation of His12 or His119, which are present in the active site of the enzyme. Nevertheless, RNase A is very active under very different conditions and difficult to inactivate.

At low salt concentrations (up to 100mM NaCl), RNase A cleaves single- and double-stranded RNA and RNA in RNA - DNA hybrids. Under high salt concentrations (>300mM NaCl) single-stranded RNA is cleaved only. To remove the enzyme from samples, it has to be digested by Proteinase K ([M3036](#), [M3037](#)) (frequently, SDS at a final concentration of 0.6% is added) and several phenol extractions are required.

Applications

- Enzymatic manipulation of DNA and RNA
- Minipreps of plasmid DNA
- In-situ hybridisation of cellular RNA
- Removal of RNA from plasmid preparations.

Specifications

Molecular weight:	approx. 13.7 kDa
CAS-Number:	[9001-99-4]
Activity:	min. 80 U/mg (Kunitz)
Extinction coefficient:	E(1%) = 7.1 (280 nm)
Isoelectric point:	pI = 9.6.
Optimal temperature:	60°C (activity range of 15°C - 70°C)
Optimal pH:	7.6 (activity range of pH 6 - 10)
Activators of RNase A:	potassium and sodium salts

Unit definition One unit is defined as the amount of enzyme necessary to hydrolyse RNA to yield a velocity constant, $k = 1$, at 25°C and pH5.0 (also known as Kunitz-Unit).

Stability / Behaviour

- RNase A aggregates while lyophilisation and storage.
- The enzyme shows a strong affinity to glass surfaces, binding strongly to glass.
- At neutral pH (for example PBS pH7.4) and at high concentrations (>10mg/mL) the enzyme will precipitate.

Storage:

-20° C.

Stored as a solution at +2° C to +8° C RNase A will be stable for several weeks.

Stored as a solution at -20° C RNase A will be stable for several years.

Lyophilised and stored at a dry place at +2° C to +8° C RNase A is stable for several years.

Stock solutions

Recommended stock solutions are from 1 - 10mg/mL in 10mM Tris/HCl, pH7.4, 15mM NaCl or in 10mM Tris/HCl, pH7.4, 1mM EDTA, pH8.0 (TE-buffer).

Working concentrations

The recommended working concentration is from 0.1 to 10µg/mL.

10µg/mL (Removal of RNA from plasmid preparations: 1 hr, RT)

100ng/mL (Preparation of "blunt ends" of double-stranded cDNA).