



# SNP Pol DNA Polymerase

## 2X PCR Master Mix

Product	Cat#	Package size
SNP Pol DNA Polymerase 2X PCR Master Mix	M3061.0100	1.25mL
SNP Pol DNA Polymerase 2X PCR Master Mix	M3061.0500	5 x 1.25mL

### Product description

The included SNP Pol DNA Polymerase is an engineered DNA polymerase variant, specially developed for all approaches in which high single nucleotide discrimination is required, for instance in allele-specific amplifications (ASA) by PCR and primer extensions or methylation-specific PCRs (MSP), HLA genotyping or the analysis of single CpG methylation sites. SNP Pol DNA Polymerase efficiently amplifies from primers that are matched at the 3'-end and discriminates primers that are mismatched! The kit comes together with an optimized buffer system. SNP Pol DNA Polymerase 2X PCR Master Mix ensures reproducible results, significantly reduces set-up times and the risk of pipetting mistakes. The SNP Pol DNA Polymerase 2X PCR Master Mix can also be used for realtime cycling when adding a suitable realtime PCR fluorescence dye.

Carry-over contamination in PCR can be controlled through an enzymatic DNA digestion by adding [Uracil-DNA glycosylase \(UDG\)](#).

**SNP Pol DNA Polymerase efficiently discriminates primers, which have a mismatch at the 3'-end.**

The engineered polymerase is also available as [SNP PolTaq DNA Polymerase](#) version which has a 5'-3'-nuclease activity and therefore can be used for hydrolysis probe-based assays (Taqman® probes, molecular Beacons, etc.).

### Contents

SNP Pol DNA Polymerase 2X PCR Master Mix contains all components necessary for PCR based genotyping including an engineered DNA polymerase, an optimized reaction buffer, and ultrapure [dNTPs \(M3015 from Genaxxon\)](#). A hot-start formulation of the included DNA polymerase prevents false amplification. Only target specific primers and template need to be added. Please note that the mix contains a balanced mix of dUTP and dTTP guaranteeing high activity and giving the possibility for UDG based carry-over contamination prevention.

### Product Specifications

Concentration:	2-time master mix
Extension rate:	0,5 to 1 kb/min. at 72 °C
5'-3' exonuclease activity:	No
3'-5' exonuclease activity:	No
Nuclease contamination:	No
Protease contamination:	No

### Unit definition

One unit of SNP Pol DNA Polymerase is defined as the amount of enzyme that incorporates 10nmol of dNTP's into acid-insoluble fraction in 30 minutes at 72 °C under standard assay conditions.

### Quality Control

SNP Pol DNA Polymerase 2X PCR Master Mix is tested for successful ASA performance detecting a genomic SNP (rs72921001) in HeLa genomic DNA. The activity of SNP Pol DNA polymerase is monitored and adjusted to a specific DNA polymerase activity using an artificial DNA template and a DNA primer.

Enzyme concentration has been determined by protein-specific staining.

Please inquire more information at [info@genaxxon.com](mailto:info@genaxxon.com) for the lot-specific concentration.

Enzyme concentration is determined by protein-specific staining.



## Application

- SNP-detection by allele-specific amplification (ASA) / Allele-specific PCR
- Methylation specific PCR (MSP)
- HLA genotyping
- Multiplex PCR

## Stability

Genaxxon bioscience SNP Pol DNA Polymerase 2X Master Mix is shipped on wet ice but retain full activity at RT (+2°C to +8°C) for at least 6 9 months.

SNP Pol and SNP PolTaq DNA Polymerase, including buffers and reagents, should be stored immediately upon receipt at -20°C. When stored under these conditions and handled correctly, these products can be kept at least until the expiration date (see tube label) without showing any reduction in performance. The Genaxxon bioscience SNP Pol DNA Polymerase 2X Master Mix can also be stored at +2°C - +8°C up to 9 months.

## Product Use Limitations

The product (SNP Pol DNA Polymerase and SNP PolTaq DNA Polymerase) are covered by a pending patent application. It is the purchaser's own internal research use and may not be resold, modified or used for production and commercial purposes of any kind.

For information on obtaining additional rights, please contact Genaxxon bioscience ([info@genaxxon.com](mailto:info@genaxxon.com)).

This product is for research use only and may be used in vitro experiments only.

## Safety information

When working with chemicals, always wear a suitable lab coat, disposable gloves, and protective goggles.

This product does not require a Material Safety Data Sheet because it does neither contain more than 1% of a component classified as dangerous or hazardous nor more than 0.1% of a component classified as carcinogenic. However, we generally recommend, when working with chemicals, always wear a suitable lab coat, disposable gloves, and protective goggles.

Genaxxon bioscience takes no liability for damage resulting from handling or contact with this product.

More information can be found in the REGULATION (EC) No. 1272/2008 OF THE EUROPEAN PARLIAMENT AND THE COUNCIL or contact Genaxxon bioscience ([info@genaxxon.com](mailto:info@genaxxon.com))

## PCR Protocol Part

### SNP Pol DNA Polymerase 2x Master Mix PCR Protocol

This protocol serves as a guideline for PCR amplification. Optimal reaction conditions such as incubation times, temperatures, and amount of template DNA may vary and must be determined individually.

#### Important notes before getting started

- Set up all reaction mixtures in an area separate from that used for DNA preparation or PCR product analysis.
- Spin down and mix all solutions carefully before use.
- SNP Pol DNA Polymerase 2X PCR Master Mix buffer is optimized for short amplicon length (about 60 - 200bp), but also longer amplicons are possible. The addition of additional MgCl<sub>2</sub> (0.5 - 1.5mM) might be needed in case of longer amplicons (>500bp).
- The SNP Pol DNA Polymerase 2X PCR Master Mix can also be used for realtime cycling when adding a suitable realtime PCR fluorescence dye).
- SNP Pol DNA polymerase is a nuclease deficient DNA polymerase and cannot be used together with probe-based assays. For probes our SNP PolTaq DNA polymerase (M3025) is necessary.
- Use disposable tips containing hydrophobic filters to minimize cross-contamination.

#### Procedure

1. **Thaw SNP Pol DNA Polymerase 2X PCR Mastermix, and 25mM MgCl<sub>2</sub>** (if required) at RT or on ice. Keep the solutions on ice after complete thawing. Mix well before use to avoid localized differences in salt concentration.
2. A negative control (without template DNA) should always be included in every experiment. Keep the master mix on ice. **Note:** The Mg<sup>2+</sup> concentration provided by the supplied PCR-buffer will produce satisfactory results in most cases. However, in some cases (amplicon >500bp), reactions may be improved by increasing the final Mg<sup>2+</sup> concentration according to Table 2.

**Table 1: PCR reaction components (25µL PCR reaction)**

Components	Components	Final concentration
SNP Pol 2X PCR Master Mix	12.5 µL	1X
primer forward (10µM)*:	0.5µL	0.2 µM (0.05 - 1.0µM)
primer reverse (10µM):	0.5µL	0.2 µM (0.05 - 1.0µM)
Template DNA	x µL	<10ng plasmid DNA or <500ng genomic DNA
sterile, bidistilled water	x µL	up to 50µL total volume
25mM MgCl <sub>2</sub> **	x µL	only if necessary

\* Primers should ideally have a GC content of 40-60%  
 \*\* should only be necessary for amplicon lengths of >500bp

**Table 2: Final MgCl<sub>2</sub> concentration in a 25µL reaction**

Additional MgCl <sub>2</sub> conc. in reaction (mM)	0.5	1.0	1.5	2.0
Additional volume of 25mM MgCl <sub>2</sub> per 50µL reaction (µL)	0.5	1	1.5	2

**Note:** The optimal Mg<sup>2+</sup> concentration should be determined empirically but in most cases the supplied 10X buffer will produce satisfactory results.

3. **Mix the master mix thoroughly**, and dispense appropriate volumes into PCR tubes. Mix gently, for example, by pipetting the master mix up and down. It is recommended that PCR tubes are kept on ice before placing in the thermal cycler.
4. **Add template DNA** (<1µg/reaction) to the individual tubes containing the master mix. For RT-PCR, add an aliquot from the reverse transcriptase reaction. The volume added should not exceed 10% of final PCR volume.
5. When using a thermal cycler with a heated lid, do not use mineral oil. Proceed directly to step 6. Otherwise, overlay with approximately 50µL - 100µL mineral oil.
6. Program the thermal cycler according to the manufacturer's instructions. A typical PCR cycling program is outlined in Table 3. For maximum yield and specificity, temperatures and cycling times should be optimized for each new target or primer pair.

**Table 3: PCR conditions (Thermal cycler)**

Step	time	temperature	comments
Initial denaturation:	2 min.	95 °C	
3-step cycling			
Denaturation:	15 sec.	95 °C	
Annealing:	10 sec.	54 - 72 °C	approximately 3 °C to 5 °C* below lower T <sub>m</sub> of primers.
Extension:	30 sec./250bp	72 °C	approx. 0.5kbp per minute extension rate
Number of Cycles	25 - 40		
Hold		<10 °C	

**Note:** After amplification, samples can be stored at +2 °C to +8 °C overnight, or -20 °C for long term storage.

**Table 4: Recommendations for Standard PCR-Primers**

<b>Length:</b>	18-30 nucleotides
<b>GC-Content:</b>	40-60%
<b>T<sub>m</sub>:</b>	Design primer pairs with similar T <sub>m</sub> values. Optimal annealing temperature may be above OR below the estimated T <sub>m</sub> . As a starting point, use an annealing temperature of 3 °C to 5 °C below T <sub>m</sub> of the primer with the lower T <sub>m</sub> -Value.
<b>Sequence:</b>	Avoid complementarities of two or more bases at the 3' ends of primer pairs. Avoid runs of 3 or more Gs or Cs at the 3' end. Avoid a 3'-end T. Avoid complementary sequences within primer and between primer pairs.

**Table 5: Migration Chart of some Gel Tracking Dyes**

Dye in agarose gel	0.5%-1.5%	2.0%-3.0%	CAS-number	Cat-No. Genaxxon
Xylene cyanol	1000bp - 4000bp	750bp - 200bp	2650-17-1	M3312
Cresol Red	2000bp - 1000bp	200bp - 125bp	62625-29-0	M3371
Bromophenol blue	500bp - 400bp	150bp - 50bp	115-39-9	M3092
Orange G	<100bp	<20bp	1936-15-8	M3180
Tartrazine	<20bp	<20bp	1934-21-0	