# Product Information

Product : Cell Culture Medium

english version : 190319

## Instructions for dissolving powder medium to prepare sterile liquid cell culture medium

### General Information

Powdered media and salts are very hygroscopic and must be stored under dry conditions. After opening the package the whole contents must be dissolved at once. If only part of the bottle is used, the residual powder has to be stored in the dark at +2°C to +8°C and under dry conditions because of its highly hygroscopic characteristics.

Reconstitute the powder to produce 1X liquid medium only, as the different amino acids may precipitate because of their low solubility. They potentially can form salts which are of low solubility in concentrated solutions. If supplements are needed, they can be added before filtration (unsterile) or after filtration (sterile).

Use bidistilled or deionized, pyrogen-free water to reconstitute powder media. For example [Cell Culture Grade Water - M6081](#) from Genaxxon.

### Procedure for preparing sterile liquid medium by sterile filtration

1. Add water to the required quantity of powdered medium (use approx. 80% of the required amount of water so as to adjust pH later). Rinse out any remaining powder from the container. Stir until completely dissolved. The temperature of the water should be between +15°C and +30°C.
2. After all powder has been dissolved, NaHCO<sub>3</sub> will be added under stirring. The value of the amount of NaHCO<sub>3</sub> necessary can be taken from the composition table of each cell culture medium.
3. Check and adjust to the desired pH (physiological optimum is pH6.8-7.2) with 1N HCl or 1N NaOH. The pH should be adjusted to about 0.2 below the wished value, as the following sterile filtration procedure will increase the pH by degassing of CO<sub>2</sub>. This precaution is not necessary if the used powder contains Hepes or other organic buffer substances.
4. After adjustment of the pH-value, water has to be added till the end volume is reached. The solution has to be mixed thoroughly and sterile filtrated immediately using a 0.22µm filter.
5. Protect the freshly prepared liquid medium from direct light and store at +2°C to +8°C.

### Procedure for preparing liquid medium from autoclavable powder medium

1. Add water to the required quantity of powdered medium (use approx. 80% of the required amount of water so as to adjust pH later). Flush out any remaining powder from the container. Stir until completely dissolved. The temperature of the water should be between +15°C and +30°C.
2. Check pH of solution. The pH should be about 4.1. If necessary adjust with 1N HCl.
3. Best prepare aliquots of 500mL. The lids of the containers should only be opened a little for autoclaving. Autoclave 15min. at 121°C. Make sure by appropriate measurement that time and temperature have been reached. Let the autoclaved medium cool down to room temperature.
4. Add a sterile 7.5% NaHCO<sub>3</sub> solution according to the list attached to this protocol. Add also other heat sensitive supplements as e.g. glutamine, serum, antibiotics, growth factors, etc. after autoclaving.
5. Check pH of solution. The physiological optimum is pH6.8-7.2. If necessary adjust pH with 1N HCl or 1N NaOH.
6. After adjusting pH add sterile water to the final volume and mix well.
7. Protect the freshly prepared liquid medium from direct light and store at +2°C to +8°C.

# Product Information

Product : Cell Culture Medium

english version : 190319

## Procedure for preparing sterile liquid insect cell culture medium by sterile filtration

1. Add **794.76mg/L Calciumchloride dihydrat e**( $\text{CaCl}_2 \times 2\text{H}_2\text{O}$ ) to the powdered medium before adding water.
2. Add water to the required quantity of powdered medium (use approx. 80% of the required amount of water so as to adjust pH later). Flush out any remaining powder from the container. Stir until completely dissolved. The temperature of the water should be between  $+15^\circ\text{C}$  and  $+30^\circ\text{C}$ .
3. After all powder has been dissolved,  $\text{NaHCO}_3$  will be added under stirring. The value of the amount of  $\text{NaHCO}_3$  necessary can be taken from the list attached to this protocol.
4. Adjust pH to 6.3 - 0.2 with 1N HCl or 1N NaOH. The pH should be adjusted to about 0.2 below the wished value, as the following sterile filtration procedure will increase pH slightly by degassing of  $\text{CO}_2$ .

The physiological pH optimum of insect cell culture medium is more acidic than the pH for mammalian cells. At pH values above approx. 6.5, an undesired precipitation will occur with insect cell culture media containing a high concentration of calcium chloride.

5. After adjustment of the pH-value, water has to be added till the end volume is reached. The solution has to be mixed thoroughly and sterile filtrated immediately.
6. Protect the freshly prepared liquid medium from direct light and store at  $+2$  to  $+8^\circ\text{C}$ .

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