

BME - media with Earle's salts

ask also for other formulations
Fragen Sie auch nach anderen
Zusammensetzungen

| | Component | mg/L | | | |
|-------------------------|---|---------|--|--------|------------|
| Anorganic salts | CaCl ₂ (anhydrous) | 200,00 | BME with Earle's salts without L-Glutamine with 2.2 g/L NaHCO ₃ | 500 ml | C4007.0500 |
| | KCl | 400,00 | | | |
| | MgSO ₄ dried | 97,70 | | | |
| | NaCl | 6800,00 | | | |
| | NaH ₂ PO ₄ x H ₂ O | 140,00 | | | |
| Other components | D(+)-Glucose anhydrous | 1000,00 | BME with Earle's salts with L-Glutamine with 2.2 g/L NaHCO ₃ | 500 ml | C4008.0500 |
| | Hepes | 0,00 | | | |
| | Phenol red | 10,00 | | | |
| Amino acids | L-Arginine x HCl | 21,00 | BME with Earle's salts with stab. L-Glutamine without NaHCO ₃ | 500 ml | C4009.0500 |
| | L-Cystine | 12,00 | | | |
| | L-Glutamine | 0,00 | | | |
| | L-Histidine base | 8,00 | BME with Earle's salts without L-Glutamine with 25 mM Hepes with 2.2 g/L NaHCO ₃ | 500 ml | C4010.0500 |
| | L-Isoleucine | 26,00 | | | |
| | L-Leucine | 26,00 | | | |
| | L-Lysine x HCl | 36,47 | | | |
| | L-Methionine | 7,50 | | | |
| | L-Phenylalanine | 16,50 | BME with Earle's salts with stab. L-Glutamine with 2.2 g/L NaHCO ₃ | 500 ml | C4273.0500 |
| | L-Threonine | 24,00 | | | |
| | L-Tryptophane | 4,00 | | | |
| | L-Tyrosine | 18,00 | BME with Earle's salts with L-Glutamine without Phenol red with 2.2 g/L NaHCO ₃ | 500 ml | C4274.0500 |
| | L-Valine | 23,50 | | | |
| D(+)-Biotin | 1,00 | | | | |
| D-Ca-Pantothenate | 1,00 | | | | |
| Cholinchloride | 1,00 | | | | |
| Folic acid | 1,00 | | | | |
| myo-Inositol | 2,00 | | | | |
| Nicotinamide | 1,00 | | | | |
| Pyridoxal x HCl | 1,00 | | | | |
| Riboflavine | 0,10 | | | | |
| Thiamine x HCl | 1,00 | | | | |