




## Genaxxon BioScience

### 5-Benzylmercaptotetrazole data sheet

**Catalog number**      **Product name**      **Structure**

S5401.0009  
S5401.0018  
S5401.0090  
S5401.0180  
S5401.0000

5-Benzylmercaptotetrazole



**Technical Data:**

MW = 192.2 g/mole  
C<sub>8</sub>H<sub>8</sub>N<sub>4</sub>S  
MP = 135 – 138°C

**Storage:**

Store at RT to avoid crystallisation of the 5-Benzylmercaptotetrazole.

**Shipment:**

RT temperature

Preparation of working solutions

Amount of Activator	Volume of dry Acetonitrile
1,75 g (9 mmole)	30 ml
3.5 g (18 mmole)	60 ml
17.3 g (90 mmole)	300 ml
34.6 g (180 mmole)	600 ml