


Genaxxon BioScience

Electrophoretic Equipment – Submarine Chambers



Modell E-H1

Modell E-H2

Modell E-H3

Cat-No.

I4001.0001

I4002.0001

I4003.0001

Abmessungen (LxWxH)

180x155x85 mm

260x215x90 mm

205x235x85 mm

Buffer volume

300 ml

800 ml

900 ml

Gel size (LxW)

83 x 70 mm

140 x 120 mm

102 x 145 mm

Gel volume

30 ml

85 ml

70 ml

Run width with one / two combs

70 / 35 mm

126 / 63 mm

90 / 45 mm

Wells per comb

6, 8, 10, 12

6, 8, 10, 12, 16, 20, 24

16, 30, 40

Voltage / Current

0-150 V / 0-100 mA

0-150 V / 0-100 mA

0-150 V / 0-100 mA

Price
221,00 Euro**307,00 Euro****307,00 Euro**
Included in delivery

1 UV-transparent tray
 2 electrical cables
 2 combs
 1 bubble level

1 UV-transparent tray
 2 electrical cables
 2 combs
 1 bubble level

1 UV-transparent tray
 2 electrical cables
 2 combs
 1 bubble level