

Genaxxon BioScience

Electrophoretic Equipment – Submarine Chambers



Modell E-H4

Cat-No.

14004.0001

Abmessungen (LxWxH)

280x340x85 mm

Buffer volume

1100 ml

Gel size (LxW)

16 x 23 mm

Gel volume

180 ml

Run width with one / four

150 / 32 mm

Wells per comb

25, 50

Voltage / Current

0-150 V / 0-100 mA

Price
459,00 Euro

Modell E-H5

14005.0001

415x295x100 mm

3000 ml

25 x 20 mm

250 ml

one / multiple
 230 / 15-170 mm

2, 24, 42, 50

0-150 V / 0-100 mA

515,00 Euro

Modell E-H6

14006.0001

150x160x52 mm

140 ml

83 x 102 mm

40 ml

One / two combs
 75 / 35 mm

4, 6, 12, 16

0-200 V / 0-50 mA

229,00 Euro
Included in delivery

 1 UV-transparent tray
 2 electrical cables
 2 combs
 1 bubble level
 4 levelling feet

 1 UV-transparent tray
 2 electrical cables
 3 combs
 1 bubble level
 4 levelling feet
 2 gel casting gates

 2 electrical cables
 2 combs
 1 bubble level
 2 gel casting gates