

# Genaxxon BioScience

## Electrophoretic Equipment – Submarine Chambers



### Modell E-H7

**Cat-No.**

I4007.0001

**Abmessungen (LxWxH)**

200x192x843 mm

**Buffer volume**

200 ml

**Gel size (LxW)**

113 x 102 mm

**Gel volume**

55 ml

**Run width with one / two combs**

105 / 52 mm

**Wells per comb**

8, 13, 16

**Voltage / Current**

0-300 V / 0-100 mA

**Price**
**167,00 Euro**

### Modell E-pbO

I4008.0001

200x220x43 mm

200 ml

113 x 102 mm

55 ml

105 / 52 mm

8, 13, 16

0-130 V / 0-70 mA

**295,00 Euro**

### Modell E-pb70

I4009.0001

78x110x18 mm

40 ml

60 x 50 mm

15 ml

**One comb**

42 mm

5, 10

0-200 V / 0-50 mA

**459,00 Euro**
**Included in delivery**

 2 electrical cables  
 1 comb  
 1 bubble level

 1 UV-transparent tray  
 2 electrical cables  
 3 combs  
 1 bubble level  
 1 power supply unit

 1 UV-transparent tray  
 2 electrical cables  
 2 combs  
 1 bubble level  
 1 Controller/Timer unit