

Fit Chart PCR Plates

	I2010/I2016 Full Skirted	I2003/I2004 Non-Skirted	I2014/I2017 Semi-Skirted
Standard PCR Instruments	96 Well	96 Well	96 Well
ABI			
Veriti 384well Block			
Veriti 0.1mL 96well Block			
Veriti 0.2mL 96well Block			☑
GeneAmp® 2400			
GeneAmp® 2700/2720		☑	☑
GeneAmp® 9600		☑	☑
GeneAmp® 9700		☑	☑
GeneAmp® 9800 FAST Block			
Biometra			
Uno	☑	☑	
Uno II		☑	
Tpersonal Cyclor			
T1 Thermocycler	☑	☑	☑
T3 Thermocycler			
TGradient	☑	☑	☑
TRobot	☑	☑	
Bio-Rad			
Genecycler			
iCycler™		☑	☑
MyCycler		☑	☑
C1000	☑	☑	☑
S1000	☑	☑	☑
Bio-Rad / MJ Research			
PTC100™ (96well block only)	☑	☑	☑
PTC200™ (DNA Engine™)	☑	☑	☑
PTC220/221 (DNA Dyad™)	☑	☑	☑
PTC225 (DNA Tetrad™)	☑	☑	☑
Corbett Research			
Palm Cyclor	☑		
Eppendorf			
Master-Cyclor®			☑
Master-Cyclor® Gradient	☑	☑	
Master-Cyclor® EP Gradient M384	☑	☑	☑

	I2010/I2016 Full Skirted	I2003/I2004 Non-Skirted	I2014/I2017 Semi-Skirted
Standard PCR Instruments	96 Well	96 Well	96 Well
Ericomp			
SingleBlock System		☺	
TwinBlock System		☺	
Deltacycler I		☺	
G-Storm			
GS1, GS2, GS3	☺		
MWG			
Primus 96	☺	☺	
Primus 384			
The Q-Lifecycler	☺	☺	☺
Peqlab			
peqSTAR96			☺
SensoQuest			
Labcycler	☺	☺	
Labcycler Gradient	☺	☺	
Stratagene			
Robocycler		☺	
Gradient Cyclor	☺	☺	
Takara			
TP240			
TP3000		☺	
Techne			
TC412/512	☺	☺	
Touchgene Gradient	☺	☺	
Flexigene	☺	☺	
Genius	☺	☺	
Thermo			
PCR Express/Px2/PxE	☺	☺	
Multiblocks System & MBS	☺	☺	
Touchdown	☺	☺	
Omnigene	☺	☺	
Omni-E	☺	☺	

Fit Chart PCR Plates

	I2010/I2016 Full Skirted	I2003/I2004 Non-Skirted	I2014/I2017 Semi-Skirted
qPCR Instruments	96 Well	96 Well	96 Well
ABI			
7000/7300		👍	👍
7500/7700/7900		👍	👍
FAST 7500/7900/7900HT			
7900HT			
StepOne Plus™		👍	👍
Bio-Rad			
iCycler iQ™		👍	👍
iQ™ 4, iQ™ 5		👍	👍
MyiQ™		👍	👍
CFX96	👍		
CFX384			
Bio-Rad-MJ Research			
Opticon™	👍		
Opticon2™	👍		
Chromo4™	👍		
Eppendorf			
Mastercycler™ ep realplex	👍	👍	👍
Roche			
LC480			
Stratagene			
Mx3000P™		👍	👍
Mx3005P™		👍	👍
Mx4000™		👍	👍
Techne			
Quantica			👍

Fit Chart PCR Plates

	I2010/I2016 Full Skirted	I2003/I2004 Non-Skirted	I2014/I2017 Semi-Skirted
Sequencers	96 Well	96 Well	96 Well
ABI			
3100 Genetic Analyser		👍	👍
3130 Genetic Analyser		👍	👍
310 Genetic Analyser		👍	👍
3500 Genetic Analyser		👍	
3500XL			👍
3500			👍
3700 DNA Analyser		👍	👍
3730/3730XL DNA Analyser		👍	👍
Amersham			
MegaBACE™ 500	👍		
MegaBACE™ 1000 mark 2	👍*		
MJ Research			
BaseStation™	👍		
Transgenomic			
Wave	👍		

* Compatible with MegaBACE™ bought after July 2000.