
fon:
+49 (0)731 - 3608 123
fax:
+49 (0)731 - 3608 962
eMail:
info@genaxxon.com
internet:
www.genaxxon.com

Certificate of Analysis

Product: ITS-Solution (100X)

Catalogue No.: C4304.0005
Produced: 21.03.2007
Shelf life: 31.03.2009 at -20°C
Lot No.: PX470806
Batchsize: 300mL

filled in: Cell Culture Plastic Bottles (10mL)

Prefilter: None
Sterile Filter: Sartolab plus, 0.22µm, Lot-No. 020671

Batch content

Insulin human	0.5mg/mL (150mg per Batch)
Transferrin human	0.5mg/mL (150mg per Batch)
Sodiumselenit	0.5µg/mL (150µg per Batch)
EBSS	8.34mg/mL (2.5g per Batch)

pH: 8.5 (acceptance criteria: 7.5-8.5)

Sterility: negative (incubation at 37°C for 10 days)
negative (incubation at 20°C for 10 days)

Dr. Norbert Tröndle
Genaxxon BioScience GmbH

21.03.2007