




5 (and 6)-carboxynaphthofluorescein (mixed and single isomers) data sheet:

Catalog
number

Spectrum

Structure

S5422.0020
S5422.0100



Technical
Data:

MW = 476.44 g/mol
C₂₉H₁₆O₇
soluble in DMSO, DMF and pH > 6
Abs.: 598 nm ($\epsilon = 49,000 \text{ M}^{-1}\text{cm}^{-1}$), pH 9.0
Em.: 658 nm

Storage:

Store at RT, protect from light, especially in solution.

Shipment:

RT temperature