




fon:
+49 (0)731 - 3608 123
fax:
+49 (0)731 - 3608 962
eMail:
info@genaxxon.com
internet:
www.genaxxon.com

pUC19 DNA - undigested



DEFINITION Cloning vector pUC19 complete sequence.

2686 bp circular DNA

ACCESSION M77789

VERSION M77789.2

```

1  gcgccaata cgcaaaccgc ctctccccgc gcgttggccg attcattaat gcagctggca
61  cgacaggttt cccgactgga aagcggggcag tgagcgcaac gcaattaatg tgagttagct
121 cactcattag gcaccccagg ctttacactt tatgcttccg gctcgtatgt tgtgtggaat
181 tgtgagcgga taacaatttc acacaggaaa cagctatgac catgattacg ccaagcttgc
241 atgcctgcag gtcgactcta gaggatcccc gggtagccgag ctcgaattca ctggccgtcg
301 ttttacaacg tcgtgactgg gaaaaccctg gcgttaccga acttaatcgc cttgcagcac
361 atcccccttt cgccagctgg cgtaatagcg aagaggcccc caccgatcgc ccttcccaac
421 agttgcgag cctgaatggc gaatggcgcc tgatgcggta ttttctcctt acgcatctgt
481 gcggtatttc acaccgcata tgggtgactc tcagtacaat ctgctctgat gccgcatagt
541 taagccagcc ccgacaccgc ccaacaccgc ctgacgcgcc ctgacgggct tgtctgctcc
601 cggcatccgc ttacagacaa gctgtgaccg tctccgggag ctgcatgtgt cagaggtttt
661 caccgtcatc accgaaacgc gcgagacgaa agggcctcgt gatacgccta tttttatagg
721 ttaatgtcat gataataatg gtttcttaga cgtcaggtgg cacttttcgg ggaaatgtgc
781 gcggaacccc tatttgttta tttttctaaa tacattcaaa tatgtatccg ctcatgagac
841 aataaccctg ataaatgctt caataatatt gaaaaaggaa gagtatgagt attcaacatt

```

```

901 tccgtgtcgc ccttattccc ttttttgccg cattttgcct tcctgttttt gctcaccag
961 aaacgctggg gaaagtaaaa gatgctgaag atcagttggg tgcacgagtg ggttacatcg
1021 aactggatct caacagcggg aagatccttg agagttttcg ccccgaagaa cgttttccaa
1081 tgatgagcac ttttaaagtt ctgctatgtg gcgcggtatt atcccgtatt gacgccgggc
1141 aagagcaact cggtcgccgc atacactatt ctcagaatga cttggttgag tactcaccag
1201 tcacagaaaa gcatcttacg gatggcatga cagtaagaga attatgcagt gctgccataa
1261 ccatgagtga taacactgcy gccaaacttac ttctgacaac gatcggagga ccgaaggagc
1321 taaccgcttt tttgcacaac atgggggatc atgtaactcg ccttgatcgt tgggaaccgg
1381 agctgaatga agccatacca aacgacgagc gtgacaccac gatgcctgta gcaatggcaa
1441 caacgttgcy caaactatta actggcgaac tacttactct agcttcccgg caacaattaa
1501 tagactggat ggagggcggat aaagttgcag gaccacttct gcgctcggcc cttccggctg
1561 gctggtttat tgctgataaa tctggagccg gtgagcgtgg gtctcgcggt atcattgcag
1621 cactggggcc agatggtaag cctcccgtc tcgtagttat ctacacgagc gggagtcagg
1681 caactatgga tgaacgaaat agacagatcg ctgagatagg tgctcactg attaagcatt
1741 ggtaactgcy agaccaagtt tactcatata tactttagat tgatttaaaa cttcattttt
1801 aatttaaaaag gatctaggtg aagatccttt ttgataatct catgaccaa atccctaac
1861 gtgagttttc gttccactga gcgtcagacc ccgtagaaaa gatcaaagga tcttcttgag
1921 atcctttttt tctgcgcgta atctgctgct tgcaaacaaa aaaaccaccg ctaccagcgg
1981 tggtttgttt gccggatcaa gagctaccaa ctctttttcc gaaggtaact ggcttcagca
2041 gagcgcagat accaaatact gttcttctag tgtagccgta gttaggccac cacttcaaga
2101 actctgtagc accgcctaca tacctcgctc tgctaactct gttaccagtg gctgctgcca
2161 gtggcgataa gtcgtgtctt accgggttg actcaagacg atagttaccg gataaggcgc
2221 agcggtcggg ctgaacgggg ggttcgtgca cacagcccag cttggagcga acgacctaca
2281 ccgaactgag atacctacag cgtgagctat gagaaagcgc cacgcttccc gaaggagaa
2341 aggcggacag gtatccggta agcggcaggg tcggaacagg agagcgcacg agggagcttc
2401 cagggggaaa cgcctggtat ctttatagtc ctgctcgggt tcgccacctc tgacttgagc
2461 gtcgattttt gtgatgctcg tcaggggggc ggagcctatg gaaaaacgcc agcaacgcgg
2521 cctttttacg gttcctggcc ttttgctggc cttttgctca catgttcttt cctgcgttat
2581 cccctgattc tgtggataac cgtattaccg cctttgagtg agctgatacc gctcgcgcga
2641 gccgaacgac cgagcgcagc gagtcaagtga gcgaggaagc ggaaga

```