

Schneider's Drosophila

ask also for other formulations
Fragen Sie auch nach anderen
Zusammensetzungen

	Supplement	mg/L			
Anorganic salts	KCl	1600,00	Schneider's Drosophila with L-Glutamine with 0.4 g/L NaHCO ₃	500 ml	C4193.0500
	KH ₂ PO ₄	450,00			
	MgSO ₄ (dried)	1810,00	Schneider's Drosophila without L-Glutamine with 0.4 g/L NaHCO ₃	500 ml	C4192.0500
	NaCl	2100,00			
	Na ₂ HPO ₄	700,00			
Other components	DL-Malic acid	600,00	Schneider's Drosophila without L-Glutamine without L-Methionine with 0.4 g/L NaHCO ₃	500 ml	C4277.0500
	Succinic acid	60,00			
	Fumaric acid	60,00	Schneider's Drosophila without L-Glutamine without L-Methionine with 0.4 g/L NaHCO ₃	500 ml	C4277.0500
	D(+)-Glucose anhydrous	2000,00			
	Yeast extract	2000,00			
	a-Ketoglutaric acid sodium salt	402,66			
	D-(+)-Trehalose x 2H ₂ O	2210,00			
Amino acids	β-Alanine	500,00	Schneider's Drosophila without L-Glutamine without L-Methionine with 0.4 g/L NaHCO ₃	500 ml	C4277.0500
	L-Arginine Base	600,00			
	L-Aspartic acid	400,00			
	L-Cysteine	60,00			
	L-Cystine	16,60			
	L-Glutamine	1800,00			
	L-Glutamic acid	800,00			
	Glycine	250,00			
	L-Histidine Base	400,00			
	L-Isoleucine	150,00			
	L-Leucine	150,00			
	L-Lysine x HCl	1650,00			
	L-Methionine	150,00			
	L-Proline	1700,00			
	L-Serine	250,00			
	L-Threonine	350,00			
	L-Tryptophane	100,00			
L-Tyrosine	500,00				
L-Valine	300,00				