

**Schneider's Drosophila**

ask also for other formulations  
 Fragen Sie auch nach anderen  
 Zusammensetzungen

	<b>Supplement</b>	mg/L			
<b>Anorganic salts</b>	KCl	1600.00	Schneider's Drosophila with L-Glutamine with 0.4 g/L NaHCO <sub>3</sub>	500mL	C4193.0500
	KH <sub>2</sub> PO <sub>4</sub>	450.00			
	MgSO <sub>4</sub> (dried)	1810.00			
	NaCl	2100.00			
	Na <sub>2</sub> HPO <sub>4</sub>	700.00			
	CaCl <sub>2</sub>	600.00	Schneider's Drosophila without L-Glutamine with 0.4 g/L NaHCO <sub>3</sub>	500mL	C4192.0500
<b>Other components</b>	DL-Malic acid	600.00			
	Succinic acid	60.00	Schneider's Drosophila without L-Glutamine <b>without L-Methionine</b> with 0.4 g/L NaHCO <sub>3</sub>	500mL	C4277.0500
	Fumaric acid	60.00			
	D(+)-Glucose anhydrous	2000.00			
	Yeast extract	2000.00			
	a-Ketoglutaric acid sodium salt	402.66			
	D-(+)-Trehalose x 2H <sub>2</sub> O	2210.00			
<b>Amino acids</b>	β-Alanine	500.00			
	L-Arginine Base	600.00			
	L-Aspartic acid	400.00			
	L-Cysteine	60.00			
	L-Cystine	16.60			
	L-Glutamine	1800.00			
	L-Glutamic acid	800.00			
	Glycine	250.00			
	L-Histidine Base	400.00			
	L-Isoleucine	150.00			
	L-Leucine	150.00			
	L-Lysine x HCl	1650.00			
	L-Methionine	150.00			
	L-Proline	1700.00			
	L-Serine	250.00			
	L-Threonine	350.00			
L-Tryptophane	100.00				
L-Tyrosine	500.00				
L-Valine	300.00				