

# Genaxxon BioScience

## Serafree 1 - Serumfree allround medium

fon:  
+49 (0)7357 - 91 63 77  
fax:  
+49 (0)7357 - 91 63 78  
eMail:  
info@genaxxon.com  
internet:  
www.genaxxon.com

### Properties and Applications

Serafree 1 is a defined serum free medium for the cultivation of a variety of cells. This so-called allround medium features Iscove's Medium as the basic medium. Transferrin, Bovine Serum Albumine (BSA), Cholesterine, Lipids and a number of trace elements are added.

