



GENAXXON

b i o s c i e n c e

Analysenzertifikat Certificate of Analysis

10-fach Tris-Borat-EDTA Fertigpuffer, pH7.4
10-time Tris-Borate-EDTA ready-to-use buffer, pH7.4

Kat. Nr.: / Cat. No.:
Lot / Contr. No.:

D2030.1010
162508

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		Weißes Pulver white powder	entspricht corresponds
Gehalt		0,89M Tris 0,89M Borsäure	Entspricht
Assay		0,02M EDTA 0.89M Tris 0.89M Boric Acid 0.02M EDTA	corresponds
Gewicht Weight	Pro Beutel Per bag	170,02g +/- 1% 170.02g +/- 1%	Entspricht corresponds
Löslichkeit: solubility:		Löslich in Wasser soluble in water	entspricht corresponds
Spezifikation: Specification:	Visuelle Beurteilung der Lösung. Visuell check of solution.	Klare Lösung clear solution	Entspricht corresponds
pH: pH:	Der 1-fachen Lösung of 1-time solution	8,3+/- 0,05 8.3 +/- 0.05	entspricht corresponds
Haltbarkeit: Stability:		Mind. 3 Jahre, falls bei RT gelagert Min. 3 years if stored at RT.	Juni 2013 June 2013

Substance shipped and stored at RT / wird bei RT verschickt und gelagert.

Safety Data: not hazardous / nicht gefährlich.

References:

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

Genaxxon bioscience GmbH

Datum/date: 03.06.2010

Dr. Norbert Tröndle

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This document has been generated by computer and does not contain a signature.

Genaxxon bioscience GmbH - Söflinger Str. 100 - D-89077 Ulm - Tel.: +731 3608 123 - Fax: +731 3608 962