

TC-100 Insect Medium

 ask also for other formulations
 Fragen Sie auch nach anderen
 Zusammensetzungen

	Supplement	mg/L			
Anorganic salts	CaCl ₂ anhydrous	980	TC-100 Insect Medium without L-Glutamine with 0.35 g/L NaHCO ₃	500 ml	C4195.0500
	KCl	2900,00			
	MgCl ₂ x 6 H ₂ O	2282,59			
	MgSO ₄ dried	1370,00			
	NaH ₂ PO ₄ x H ₂ O	970,00			
Other components	D(+)-Glucose anhydrous	1000,00	TC-100 Insect Medium with L-Glutamine with 0.35 g/L NaHCO ₃	500 ml	C4194.0500
	Bacto - Tryptose	2600,00			
Amino acids	L-Alanine	225,00	TC-100 Insect Medium without L-Glutamine without sulfate with 0.35 g/L NaHCO ₃	500 ml	C4292.0500
	L-Arginine Base	550,00			
	L-Aspartic acid	350,00			
	L-Asparagine x H ₂ O	391,97	TC-100 Insect Medium without L-Glutamine without L-Methionine with 0.35 g/L NaHCO ₃	500 ml	C4296.0500
	L-Cystine	20,00			
	L-Glutamine	600,00			
	L-Glutamic acid	600,00			
	Glycine	650,00			
	L-Histidine x HCl x H ₂ O	3400,00			
	L-Isoleucine	50,00			
	L-Leucine	75,00			
	L-Lysine x HCl	630,00			
	L-Methionine	50,00			
	L-Phenylalanine	150,00			
	L-Proline	350,00			
	L-Serine	550,00			
	L-Threonine	180,00			
	L-Tryptophane	100,00			
L-Tyrosine	55,00				
L-Valine	100,00				
Vitamines	p-Aminobenzoic acid	0,02			
	D(+)-Biotine	0,01			
	D-Ca-pantothenate	0,11			
	Folic acid	0,02			
	myo-Inositol	0,02			
	Nicotinic acid	0,02			
	Pyridoxin x HCl	0,02			
	Riboflavine	0,02			
	Thiamine x HCl	0,02			
	Vitamine B12	0,01			