



Certificate of Analysis

Product: Sterile filtered Water for Cell Culture in 20 L Bag

Catalogue No.: M6081.1002
Produced: 20.12.2007
Shelf life: 31.03.2009 at room temperature
Lot No.: PX20122008
No of bags: 30 x 20 L

filled in: Cell Culture Plastic Bags from Stedim (Flexboy 20 L)

Prefilter: None
Sterile Filter: Millidisk 20S, 0,22um

Conductivity: 0,050 uS/cm (< 0,06 uS/cm)

Endotoxin: < 0,005 EE/ml (< 0,01 EE/ml)

Cytotoxicity: negative (negative)

Sterility: negative (incubation at 37°C for 10 days)
negative (incubation at 25°C for 10 days)

Dr. Norbert Tröndle
Genaxxon BioScience GmbH

20.12.2007