

Fit-Chart: Standard PCR-Plates

qPCR Cyclers

Supplier	Cycler	Prod.-No. Genaxxon
ABI/ LIFE TECHNOLOGIES	7000, 7300, 7500, 7700, 7900, QuantStudio™, ViiA7™	I2014
BIO-RAD	CFX96, iCycler™ IQ, IQ™ 4, IQ™ 5, MyiQ™	I2010, I2016, I2004
BIO-RAD MJ RESEARCH	Chromo4™, Opticon™, Opticon2™	I2010, I2016, I2004
EPPENDORF	Mastercycler™ ep realplex	I2010, I2016, I2004, I2014
ROCHE	LC96	I2249
STRATAGENE/ AGILENT	Mx3000P™, Mx3005P™, Mx4000™	I2010, I2016, I2004, I2014
TECHNE	Quantica	I2014

Standard Thermal Cyclers

ABI/ LIFE TECHNOLOGIES	GeneAmp® 2700, 2720, 9600, 9700, Verity 0.2mL (96-well block)	I2004, I2014
BIOMETRA	T1 Thermocycler, TGradient, TRobot, Uno, Unoll	I2010, I2016, I2004, I2014
BIO-RAD	C1000, iCycler™, MyCycler™, S1000	I2010, I2016, I2004, I2014
BIO-RAD MJ RESEARCH	Mini Gradient, PTC100™ (96-well block only), DNA Engine™, DNA Dyad™, DNA Tetrad™	I2010, I2016, I2004, I2014
CORBETT RESEARCH	(Qiagen) Palm Cycler	I2010, I2016
EPPENDORF	MasterCycler®, MasterCycler® EP Gradient, MasterCycler® Gradient	I2010, I2016, I2004, I2014
MWG	Primus 96, The Q Lifecycler	I2010, I2016, I2004, I2014
PEQLAB	peqSTAR 96	I2010, I2016, I2004, I2014
STRATAGENE	Gradient Cycler	I2010, I2016, I2004
TECHNE	Cyclogene, Flexigene, Genius, Genius Quad, TC412/512, Touchgene, Touchgene Gradient	I2010, I2016, I2004
THERMO FISHER SCIENTIFIC	MultiBlock System, MBS, PCR Sprint	I2010, I2016, I2004

DNA Sequencers

ABI/ LIFE TECHNOLOGIES	310, 3100 Genetic Analyser, 3130 Genetic Analyser, 3500 Genetic Analyser, 3500XL, 3700, 3730, 3730XL DNA Analyser	I2004, I2014
AMERSHAM	MegaBACE™ 1000 mark 2, MegaBACE™ 500	I2010, I2016
MJ RESEARCH	BaseStation™	I2010, I2016
TRANSGENOMIC	Wave	I2010, I2016, I2004

Plate information is for guidance only; please contact Genaxxon bioscience which plate is suitable for your cycler.

Fax: +49 731 3608 962

Tel.: +49 731 3608 123

September 2014