

# Available Standard Formats for qPCR lyophilized LyoBalls

The Genaxxon LyoBalls are supplied packed in 2 standard formats (see below).

- Please choose your preferred format and version.
- Other plate and tube versions on request.

## Standard format 1 (PCR plates with 0.2mL or 0.1mL wells):

96 LyoBalls pre-dispensed in a Break-a-way PCR Plate with standard wells (0.2mL) or as low profile plates (0.1mL) sealed with 8 cap strips.

The standard colour of the plate frame is blue. The wells are clear, transparent polypropylene.

Figure 1: Break-a-way plate



Figure 2: Easy to divide Break-a-way plate




Figure 3: 8 cap strip for sealing the Break-a-way plate



## Standard format 2 (PCR strips with 0.2mL or 0.1mL tubes):

96 LyoBalls pre-dispensed in 12 x 8 PCR tube strips with standard wells (0.2mL) or as 8 PCR tube strips with low profile tubes (0.1mL). Each version with individually attached caps.

The tubes are from clear, transparent polypropylene.

Figure 4: 8 tube PCR strips with individually attached caps.

