

ROX level and applied real-time instrument

	Bio-Rad					Roche				Qiagen/Corbett		Life Technologies					Agilent Technologies			Thermo	Cepheid	Bio Molecular Systems	
	CFX96Touch™ & CFX384 Touch™	CFX Connect™	Opticon® 2	Chromo4™	iCycler iQ™ & MyiQ™	Lightcycler® 480	Lightcycler® 1536	Lightcycler® Nano	Lightcycler® 96	Rotor-Gene Q	Rotor-Gene 6000	7500, 7500 Fast	Vii™A7	QuantStudio™ 12K Flex	7000 7300, 7700, 7900, 7900HT	StepOne™, StepOnePlus™	Mx3000P™, Mx3005P™	Mx4000™	AriaMx	PikoReal™	SmartCycler®	Mic qPCR cyclor	
RealQ Plus Green, High ROX™															x	x							
RealQ Plus Green, Low ROX™											x	x	x				x	x	x				
RealQ Plus Green, Without ROX™	x	x	x	x	x	x	x	x	x	x	x			x*						x	x	x	
RealQ Plus Probe, High ROX™															x	x							
RealQ Plus Probe, Low ROX™											x	x	x				x	x	x				
RealQ Plus Probe, Without ROX™	x	x	x	x	x	x	x	x	x	x	x			x*						x	x	x	

*For OpenArray® experiments