

RealTime Mastermix Selection Chart

	BioRad				Roche				Qiagen/Corbett		Life Technologies				Agilent Technologies			Thermo	Cepheid			
	CFX96 Touch™ & CFX384 Touch™	CFX Connect™	Option® 2	Chromoe4™	iCycler iQ™ & MyiQ™	Lightcycler® 480	Lightcycler® 1536	Lightcycler® Nano	Lightcycler® 96	Rotor-Gene Q	Rotoer-Gene 6000	7500, 7500 Fast	Vii™A7	QuantStudio™ 12K Flex	7000, 7300, 7700, 7900, 7900HT	StepOne™, StepOnePlus™	Mx3000™	Mx3005P™	Mx4000™	PikoReal™	SmartCycler®	
GreenMastermix, High ROX™ (M3052)															X	X						
GreenMastermix, Low ROX™ (M3011)												X	X	X			X	X	X			
GreenMastermix, No ROX™ (M3023)	X	X	X	X	X	X	X	X	X	X	X			X*						X	X	
ProbeMastermix, High ROX™ (M3010)															X	X						
ProbeMastermix, Low ROX™												X	X	X			X	X	X			
ProbeMastermix, No ROX™ (M3045)	X	X	X	X	X	X	X	X	X	X	X			X*						X	X	

* For openArray® experiments