

## SDS Solution (5%)

|                                | Cat# | M3315.0250 | M3315.1000 |
|--------------------------------|------|------------|------------|
| Component                      |      |            |            |
| 5% SDS dissolved in pure water |      | 250mL      | 1L         |

### Description

The 5% SDS solution is an aqueous solution containing 5% detergent to clean optical lenses of machines and instruments. After using our 5% SDS solution the lenses have to be rinsed with our cleansing solution M6080.

### Suitability

For research usage only.

Not approved for human or animal diagnostic or therapeutic procedures. Not for consumption.

|           |                             |
|-----------|-----------------------------|
| Content:  | 50g SDS/1L water            |
| Purity:   | min. 99% of used substances |
| pH value: | 7.5 to 9.0                  |
| Storage:  | store at RT after receipt   |