

# Realtime PCR/qPCR

**Embrace Precision in Your qPCR!**

**Your Research, Our Commitment!**

**Our expertly designed qPCR mastermixes deliver unparalleled accuracy and consistency.**

**Achieve reliable results for your research every time.**

**We Aim for Your Success.**

**[www.genaxxon.com](http://www.genaxxon.com)**

# Our Best Selling qPCR Mastermixes, Customized for Your Needs

Choose **GreenMasterMix** or **ProbeMasterMix**  
with or without an intercalating fluorescent dye,  
and select options with or without ROX.



- GREENMASTERMIX NO ROX (M3023)
- GREENMASTERMIX LOW ROX (M3011)
- GREENMASTERMIX HIGH ROX (M3052)
  
- PROBEMASTERMIX NO ROX (M3045)
- PROBEMASTERMIX LOW ROX (M3031)
- PROBEMASTERMIX HIGH ROX (M3010)

Our optimized mastermixes are ideal for realtime PCR in block systems, delivering all necessary components for flawless qPCR. Simply add your primer and template DNA.

## Also available as FAST-qPCR master mixes

- GREENMASTERMIX FAST NO ROX (M3150)
- GREENMASTERMIX FAST LOW ROX (M3154)
- GREENMASTERMIX FAST HIGH ROX (M3228)
  
- PROBEMASTERMIX FAST NO ROX (M3144)
- PROBEMASTERMIX FAST LOW ROX (M3152)
- PROBEMASTERMIX FAST HIGH ROX (M3155)



INTERESTED?



Shop  
Now

[www.genaxxon.com](http://www.genaxxon.com)

# From Any Platform to Any Sample – Reliable qPCR Solutions


## ANY PLATFORM

### GNX VERSA qPCR MASTERMIX 2X (M3212)

#### One qPCR Master Mix for All Platforms!

Forget the hassle of selecting different qPCR master mixes for various platforms.

With the GNX Versa qPCR Master Mix 2X, you only need one product that's compatible with most real-time qPCR instruments. Designed for streamlined setup and consistent results, it's the ideal choice for labs looking to simplify their workflows without compromising on performance.



## ANY SAMPLE

### GNX FOOD PROBE qPCR MASTERMIX 2X (M3048)

#### One qPCR Master Mix for the Toughest Samples!

Simplify your workflow without compromising precision. The Genaxxon GNX Food Probe qPCR MasterMix 2X is expertly designed to tackle even the most challenging samples, such as environmental or food matrices with high levels of inhibitors.

INTERESTED?



shop  
Now

We Aim for Your Success.


[www.genaxxon.com](http://www.genaxxon.com)

# Your Multiplexing:

## Mastering Precision!

### MULTIPLEX HS MASTERMIX 2X (M3013)

Optimized ready-to-use mixture for probe-based assays such as TaqMan®, Beacons and MGBs. It contains an antibody-inhibited fast HotStart Taq DNA Polymerase, dNTPs and MgCl<sub>2</sub> combined in an optimized buffer system for qPCR applications. Just add primers, probe and template DNA / cDNA.




M 10-plex Single PCR bands 10-plex M  
10-plex Multiplex PCR with Genaxxon's Multiplex HS Mastermix 2X

### 5X qPCR MULTIPLEX MASTERMIX (M3024)

For robust multiplexing qPCR with all components for rapid, sensitive and reproducible results. The optimized DNA polymerase with hot-start function and an optimized buffer including our ultrapure dNTPs are the key components of the mastermix.

Applicable for a wide range of templates.



INTERESTED?



Shop  
Now

**GENAXXON**  
bioscience

Söflinger Str. 70  
89077 Ulm

Tel +49(0)731 - 3608-123

Fax +49(0)731 - 3608-962

E-Mail [info@genaxxon.com](mailto:info@genaxxon.com)

printed on recycled paper