



## Certificate of Analysis

fon:  
+49 (0)731 - 3608 123  
fax:  
+49 (0)731 - 3608 962  
eMail:  
info@genaxxon.com  
internet:  
www.genaxxon.com

**Product:** PLP (178 - 191)  
**Sequence:** NWTTCQSIAPFSK  
**Synonyms:** Asn-Thr-Trp-Thr-Thr-Cys-Gln-Ser-Ile-Ala-Phe-Pro-Ser-Lys  
Myelin proteolipid protein (PLP) residues 139-151  
**Catalogue No.:** P2274.9501  
**Produced/Re-Test:** 11.03.2013  
**Shelf life:** March 2015 at 4°C  
**Lot No.:** J14b  
**Appearance:** lyophilized trifluoroacetate salt (crystalline white powder)  
**delivered:** filled in glass vial.  
**Purified:** by RP-HPLC  
**Purity:** >95% (HPLC)  
**MW:** 1582.8g/mol / C<sub>70</sub>H<sub>105</sub>N<sub>18</sub>O<sub>22</sub>S  
**CAS-No.:** not available  
**Reconstitute:** 2mg/mL in water.

Substance shipped at RT, please store at 4°C / wird bei RT verschickt, bitte bei 4°C lagern.

**Safety Data:** not hazardous / nicht gefährlich



Dr. Norbert Tröndle  
Genaxxon bioscience GmbH

11.03.2013

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.  
This document has been generated by computer and does not contain a signature.