



Certificate of Analysis

fon:
+49 (0)731 - 3608 123
fax:
+49 (0)731 - 3608 962
eMail:
info@genaxxon.com
internet:
www.genaxxon.com

Product: OVA (257-264)
Sequence: SIINFEKL
Synonyms: Ser-Ile-Ile-Asn-Phe-Glu-Lys-Leu, Ovalbumin (257 - 264) chicken, Ova 257-264, class I (Kb)-restricted peptide epitope of OVA
Catalogue No.: P2275.7005
Produced/Re-Test: 11.01.2015
Shelf life: Jan. 2018 at 4°C
Lot No.: G13b
Appearance: lyophilized trifluoroacetate salt (crystalline white powder)
delivered: filled in glass vial.
Purified: by RP-HPLC
Purity: >70% (HPLC)
MW: 963.13g/mol / C₄₅H₇₄N₁₀O₁₃
CAS-No.: 138831-86-4
Reconstitute: 2mg/mL in water.

Substance shipped at RT, please store at 4°C / wird bei RT verschickt, bitte bei 4°C lagern.

Safety Data: not hazardous / nicht gefährlich



MW (calc.) 963
[M+H]⁺ 964
[M+2H]²⁺ 483

Dr. Norbert Tröndle
Genaxxon bioscience GmbH

11.01.2015

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This document has been generated by computer and does not contain a signature.