

Genaxxon bioscience GmbH • Söflinger Str. 70 • D – 89077 Ulm

Summary of the analytical results of the peptide synthesis

Stock product P2277

Sequenz	Name	Batch-No.:	Appearance	MW (calc.) g/mol	Purity % (HPLC (214nm))	Next retest
H-GILGFVFTL-OH	Influenza A MP (58-66) HLA-A*02:01	0205V07	White powder	966.4	97.2%	05/2028

Counterion: TFA (Trifluoroacetate)

HPLC-Chromatogram(s) and MS spectra are part of the shipment.

These products are for research and development usage only!

Please contact Genaxxon bioscience (info@genaxxon.com) for information if you plan to use the products for other purposes.

Note: Please store peptides at -20°C in a dry location.

Note: Do not open the tube/bottle before peptides have been warmed up to ambient temperature.

Complaints can be generally only accepted within the first week after receipt of the peptides.

Complaints due to improper handling or storage cannot generally be accepted.



Zertifikat-Registrier-Nr.
12 100 39900 TMS

GENAXXON bioscience GmbH
 Söflinger Str. 70 • 89077 Ulm • Geschäftsführer: Dr. Norbert Tröndle
 Amtsgericht: 89014 Ulm • HRB 641623
 U.St.-Id.Nr. DE 220 603 767 • Deutsche St.-Nr. 88002/33658

We Aim for Your Success.

www.genaxxon.com