

Summary of the analytical results of the peptide synthesis

Peptide: P2286 - Nona-Arginine - Arg9

Cat no.	Sequenz	Name	N-Term	C-Term	Batch No.:	MW (calc.) g/mol	Purity % (HPLC (214nm))
P2286	RRRRRRRRR	Nona-Arginine - Arg9	H	OH	1006X01	1423.8	96.7

Counterion: TFA (Trifluoroacetate)

HPLC-Chromatogram(s) and MS spectra are part of the shipment.

These products are for research and development usage only!

Please contact Genaxxon bioscience (info@genaxxon.com) for information if you plan to use the products for other purposes.

Note: Please store peptides at -20°C in a dry location.

Note: Do not open the tube/bottle before peptides have been warmed up to ambient temperature.

Complaints can be generally only accepted within the first week after receipt of the peptides.

Complaints due to improper handling or storage cannot generally be accepted.

