

Genaxxon bioscience GmbH • Söflinger Str. 70 • D – 89077 Ulm

- Summary of the analytical results of the peptide synthesis - for stock product MAGE-A4 (230-239) peptide - >95%

Sequenz	Name	Cat No.:	Batch No.:	MW (theo.) g/mol	Purity % (HPLC)
GVYDGREHTV	MAGE-A4 (230-239)	P2969	3003U03	1132.2	99.2

Counterion: TFA (Trifluoracetate)

HPLC-Chromatogram(s) and MS spectra are part of the shipment.

These products are for research and development usage only!

 Please contact Genaxxon bioscience (info@genaxxon.com) for information if you plan to use the products for different purposes.

Note: Please store peptides at a dry place at -20 °C.

Note: Please open the tube/bottle not before peptides have been warmed up to ambient temperature.

Complaints can be generally only accepted within the first week after receipt of the peptides.

Complaints due to improper handling or storage can generally not be accepted.

