



Analysezertifikat
Certificate of Analysis

ExoScience EVs – lyophilisierte, extrazelluläre Vesikel aus menschlichen MSC	Kat. Nr.: / Cat. No.:	C4163
ExoScience EVs - lyophilized extracellular vesicles from human MSC	Lot / Contr. No.:	260618001

lyophilized extracellular vesicles from human MSC

Parameter Parameter	Methode Method	Spezifikation Specification	Ergebnis Result
Aussehen: Appearance:		Lyophilisiertes, weißes Pulver Lyophilized, white powder	entspricht corresponds
Endotoxin:	Kinetic Turbidimetry / LAL Assay	<5 EU/mL	entspricht corresponds
Mycoplasmen:: Mycoplasma:	qPCR	ND	ND
Gesamtprotein: Total protein:	Bradford Assay	>5mg/mL	<6mg/mL
Sterilität: Sterility:	Membranfiltration Membrane filtration	ND	ND
Nächste QC: Next retest:		06/2027	
Partikel pro mL Particle per mL	Nanoparticle Tracking Analysis (NTA)	>0,5x 10E10 Partikel/mL >0.5x 10E10 particles/mL	2,05x 10E10 +/-0,34x 10E10 2.05x 10E10 +/-0.34x 10E10



Versand: nicht gekühlt. / Lagerung des lyophilisierten Produktes bei +2°C bis +25°C
 Shipment: not cooled. / Storage of the lyophilized product at +2°C to +25°C.

Storage after reconstitution in deionized water

After reconstitution of the lyophilized extracellular vesicles, we recommend storing at +2°C to +8°C for only one week!

For long-term storage we recommend preparing different aliquots and store them at -20°C. Just thaw the tube you are going to use each time. We also suggest not to store for too long at all as it is not predictable when the proteins inside exosomes will degrade.

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This document has been generated by computer and does not contain a signature.

Analysezertifikat **Certificate of Analysis**

ExoScience EVs – lyophilisierte, extrazelluläre Vesikel aus menschlichen MSC	Kat. Nr.: / Cat. No.:	C4163
ExoScience EVs - lyophilized extracellular vesicles from human MSC	Lot / Contr. No.:	260618001

lyophilized extracellular vesicles from human MSC

Safety Data: not hazardous

Usage

For research usage only!

Not for human or animal diagnostic or therapeutic usage/applications. Not for consumption.

Verwendung

Nur für Forschungszwecke!

Nicht für human- oder tierdiagnostische oder therapeutische Zwecke/Anwendungen. Nicht zum Verzehr zugelassen.

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

Genaxxon bioscience GmbH

Datum/date: 18.06.2026

Dr. Norbert Tröndle

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This document has been generated by computer and does not contain a signature.